# **EXTERNAL CARDIAC STIMULATOR DUAL-CHAMBER 1797 INFORMATION FOR SAFETY AND PERFORMANCE (EN)**

## **1 DEVICE DESCRIPTION**

FIAB 1797 external pacemaker is a compact, temporary dual chamber cardiac stimulator to be used in case of pathological alterations of the heart pulse spontaneous generation or conduction. Powered by an internal battery, it is connected to electro-catheters for temporary atrial, ventricular and AV-sequential stimulation in clinical environment.

Various stimulations modes are available, along with the high rate atrial BURST and RAMP functions.

The stimulation mode can be selected via a quick menu together with pulse rate and amplitude, input signal sensitivity (sensing threshold), PVARP, and AV-delay (AVD); the high rate stimulation modes BURST and RAMP can be selected, as well as the BURST rate and the maximal RAMP rate.

The FIAB 1797 principal characteristics are:

- Stand by, pre-stimulation modality for parameters and modes selection
- Synchronous and asynchronous stimulation
- System driven by two microprocessors, allowing high precision in selecting the stimulation parameters and a complete and continuous control of the device functioning
- Low consumption liquid crystals displays
- EMERGENCY function for a quick launch of asynchronous stimulation with safety parameters
- BURST and decreasing RAMP modes for high frequency stimulation in the treatment of tachycardia
- Backup batteries allowing main battery replacement without interrupting stimulation
- Easy and safe connection to temporary electro-catheters for transvenous epicardial or myocardial stimulation, either direct or by means of a separate patient cable or with an adaptor, when needed.

### **1.1 Safety Function**

- Self-test function
- On-line check of batteries charge level, signaling exhaustion both of the main and of the backup battery
- continuing stimulation during battery replacement
- keyboard lock preventing accidental parameters modification
- two-step switching off preventing accidental deactivation of the device
- anti-runaway protection preventing stimulation rate to go out of control
- protection against defibrillator discharge
- protection against electrostatic discharge
- protection against electromagnetic interferences
- continuous check of electrical connections, of load impedance and of stimulating current
- potentially dangerous situations trigger alarms to the user by means of acoustic signals and corresponding warnings on the display

### **1.2 Main Operation features**

- Single chamber modes: atrial AAI, AAT, A00; ventricular VVI, VVT, V00
- Dual chamber modes: DDD, DDI, D00, VDD, DAI, DVI, VAT
- Pacing and sensing LEDs, showing that electric signals are being sent (pace) or received (sense) by the device
- Low battery messages warning that the main battery has to be replaced
- Low battery messages warning that the back-up battery has to be replaced
- Upper display dedicated to setting stimulation parameters: amplitude and rate visualized and easily selectable turning the corresponding knobs
- Quick menu for selecting or changing the stimulation mode

- Lower display dedicated to supplementary menus and messages
- Screenshots of the lower display to regulate additional parameters like sensitivity and the rate dependent parameters, PVARP and AVD
- Manual or automatic selection of PVARP and AVD
- Alarm/alert messages visualized on the lower display
- LOCK/UNLOCK keyboard button, preventing accidental modifications of the parameters
- PAUSE button suspending stimulation and visualizing the patient intrinsic heart rhythm
- EMERGENCY button shifting to asynchronous dual chamber stimulation (V00) in emergencies.

### 1.3 Indications

Temporary pacemaker is used in a hospital setting for temporary single or dual chamber cardiac pacing in case of pathological alterations of the heart pulse spontaneous generation or conduction:

- cardiac arrest or complete heart block
- sinus bradycardia

### 1.4 Intended users

The 1797 dual-chamber stimulator is intended to be used by highly qualified medical personnel: specialists in cardiology or cardiac surgery, experts in cardiac stimulation. Before starting to use the device, train medical personnel on the functionality and use of the temporary pacemaker, based on needs and hospital procedures.

### 1.5 Environment of use

The DM must be used in hospital environments: cardiology, emergency room and cardiac surgery. The temporary pacemaker must be used in an environment where continuous monitoring of the patient is guaranteed, in order to ensure optimal functioning of the device and guaranteeing that the patient receives the most adequate treatment.

### 1.6 Essential Performance

The essential performances of the MD can be summarized in two fundamental functions.

- Pacing: the MD must stimulate the heart tissue (with amplitude and frequency that can be set by the user and in the chosen mode - synchronous or asynchronous)
- Sensing: the MD must recognize spontaneous cardiac activity and synchronize with it if the device has been set with ON DEMAND functionality (synchronous).

### 1.7 Clinical Benefit

Temporary cardiac stimulation is used in the treatment of heart disorders to:

- Heartbeat recovery, keeping the patient alive until spontaneous recovery or until the implantation of a permanent pacemaker
- Maintenance of the correct heart rate in case of pathological alterations of the spontaneous generation of the cardiac impulse, or stabilization of the patient from a hemodynamic point of view, i.e. guaranteeing the correct flow of blood into the patient

### 1.8 Intended purpose

A temporary pacemaker is used in connection to a system of electro-catheters, so to provide a temporary single- or dual-chamber stimulation in a hospital environment by trained personnel. The temporary pacemaker can provide either synchronous (on demand) or asynchronous, short term stimulation for therapeutic, prophylactic or diagnostic purposes, in particular:

- support, management, assessment of a patient with cardiac complications during invasive procedures
- patient support after heart surgery

## 1.9 Intended patient population

Patients who need short-term, synchronous or asynchronous cardiac electro-stimulation for therapeutic, prophylactic or diagnostic purposes.

## 2 BATTERY

The device is powered by an internal battery. Only 6LR61 type 9V batteries can be used (classification according to the international standard IEC 60086). Batteries with different dimensions or rechargeable batteries may alter the device functioning, with stimulation failure or damage to the machine. Duration of rechargeable batteries is not foreseeable and their use can cause a sudden stop of stimulation. Before use, particularly if a long time used is programmed, it is strongly advisable to replace the battery. Replacement is mandatory when the corresponding alarm is displayed on the device.

Use only fully charged new batteries, before expiration date. When the device is stored for a long time (more than 4 weeks) it is necessary to remove the battery from its housing.

## 3 POTENTIAL COMPLICATIONS – RESIDUAL RISKS

About the procedure.

Some possible complications inherent to the temporary external cardiac stimulation procedure :

- asystole after stimulation sudden interruption (e.g. due to accidental disconnection of the patient cable or of an electro-catheter, or because of a wrong setting of the device);
- signal inhibition or decay in the presence of high intensity electromagnetic interferences;
- onset of a tachyarrhythmia or acceleration of an already present one;
- dependence on the pacemaker.

Some possible complications inherent to the device itself:

- muscular stimulation causing contractions and pain;
- stimulation irregular or absent due to electro-catheter breakage or dislocation;
- loss of the signal capture due to a wrong regulation of the stimulation parameters (ineffective stimulation);
- tachycardia induced by a wrong regulation of the stimulation rate (inappropriate stimulation);
- loss of synchronism or stimulation inhibition due to a wrong regulation of the sensing parameters (Undersensing or Oversensing).

## 4 CONTRAINDICATIONS

No contraindications are known against the use of temporary cardiac stimulation as a tool for controlling and regulating the heartbeat. Nevertheless, age and clinical conditions of the patient may limit the choice of the operation mode, of the stimulation parameters and of the electro-catheters system to be employed.

Here is a list of some general caveat.

### Atrial sensing

Stimulation modes in which atrial sensing triggers a ventricular response are not indicated in the presence of rapid atrial arrhythmias like atrial fibrillation or flutter.

### Atrial stimulation

Atrial stimulation is ineffective in the presence of atrial fibrillation or flutter. Single-chamber atrial stimulation is contraindicated in the presence of A-V conduction disorders.

### Asynchronous stimulation

Asynchronous stimulation is contraindicated in the presence of intrinsic heart rhythm.

### High rate BURST and RAMP stimulation

These asynchronous functions are indicated for the treatment of atrial synchronized tachyarrhythmia.

High rate BURST stimulation is a kind of therapy to be used exclusively in atrium. If applied to ventricle it may cause lethal arrhythmias. The device is designed to apply high rate stimulation to the atrium only, but

the user has to be sure that leads are connected properly.

High rate atrial stimulation can produce tachycardia, accelerate an existing tachycardia or induce fibrillation. It has to be performed in a clinical, strictly controlled environment with continuous monitoring of ECG, blood pressure, and with an available defibrillator possibly connected to adhesive plates previously applied to the patient.

### **Competitive stimulation**

Temporary cardiac stimulation is contraindicated in the presence of another stimulation system. Do not use the device to stimulate a patient already connected to a similar system. Otherwise, the two system may enter a competition with serious consequences. For instance, a stimulation pulse can be delivered during a T-wave or a pacemaker-mediated tachycardia may arise.

## **5 WARNINGS AND PRECAUTIONS**

Necessary precautions for the correct use of the device are listed here. More indications can be found in the remaining part of these IFU.

### **WARNING!**

The 1797 device is intended to be used by highly qualified medical personnel: specialists in cardiology, cardio-surgeons, and experts in cardiac stimulation. The device has to be used in a hospital environment: cardiology and cardio-surgery departments, emergency rooms

### **WARNING!**

It is necessary to control the stimulation efficacy and the appropriate parameters setting by means of an ECG monitor. ECG patient monitoring is indispensable during the use of the device. In case of extended use of the device, it is necessary that sensing and pacing thresholds are periodically checked by a physician. A ready-to-use defibrillator must be available for emergency during the following procedures: measuring sensing and pacing thresholds; connecting the device to the electro-catheters and performing the relevant regulations; applying the atrial BURST or RAMP high rate stimulations.

### **WARNING!**

The correct determination of the sensing threshold is fundamental for the synchronous modes (AAI, VVI, DDI) to work properly. If a too high threshold is chosen, there can be a loss of capture with a consequent unnecessary, potentially dangerous stimulation (this is commonly called 'Undersensing'). A too low threshold is not recommended either, since irrelevant activities (like e.g. muscle tremors) or slight disturbances could be erroneously interpreted as spontaneous cardiac activity, inducing continuous inhibitions of synchronous stimulation; moreover, overlapped disturbances could hide the physiological signal making it interpreted as a noise with the consequence that the device will ignore it (these are the two aspects of the so-called 'Oversensing').

### **WARNING!**

Asynchronous stimulation modes (A00, V00, D00) may interfere with the heart spontaneous activity and therefore they are generally not recommended. They are only recommended in the absence of any spontaneous activity.

### **WARNING!**

Stimulation modes in which atrial sensing triggers a ventricular response are not indicated in the presence of rapid atrial arrhythmias like atrial fibrillation or flutter.

When setting a mode with atrial sensing it is necessary to assess the sensing threshold so to guarantee a sufficient safety margin. Locate the electro-catheter on the right atrium wall, oriented along the direction of the cardiac fibers and with the contacts at a distance of about 1cm from one another. It is important that the obtained sensing threshold is of at least 1mV. Set the atrial sensing threshold to a half of the measured threshold. Failing to follow this procedure may lead to the delivery of asynchronous stimulations.

A larger safety margin can be obtained by further lowering the sensing threshold. However, remember that choosing a too low value may cause a sensing of far field signals (like the R or T

waves being sensed in the atrial channel or the P wave in the ventricular channel), causing inappropriate stimulation inhibition.

**WARNING!**

Be sure that the safety margin for both atrial and ventricular sensing is adequate, otherwise the delivered stimulation could not be appropriate.

**WARNING!**

Do not use the temporary pacemaker to determine the sensing thresholds of a permanently implanted system. During the implantation of a permanent pacemaker FIAB suggests using a threshold analyzer (PSA e.g.: FIAB 2307).

**WARNING!**

Atrial stimulation is ineffective in the presence of atrial fibrillation or flutter. Single-chamber atrial stimulation is contraindicated in the presence of A-V conduction disorders.

**WARNING!**

These asynchronous functions are indicated for the treatment of atrial synchronized tachyarrhythmia.

High rate BURST stimulation is a kind of therapy to be used exclusively in atrium. If applied to ventricle it may cause lethal arrhythmias. The device is designed to apply high rate stimulation to the atrium only but the user has to be sure that leads are connected properly.

High rate atrial stimulation can produce tachycardia, accelerate an existing tachycardia or induce fibrillation. It has to be performed in a clinical, strictly controlled environment with continuous monitoring of ECG, blood pressure, and with an available defibrillator possibly connected to adhesive plates previously applied to the patient.

**WARNING!**

Temporary cardiac stimulation is contraindicated in the presence of another stimulation system. Do not use the device to stimulate a patient already connected to a similar system. Otherwise, the two system may enter a competition with serious consequences. For instance, a stimulation pulse can be delivered during a T-wave, or a pacemaker mediated tachycardia may arise.

**WARNING!**

These tests are mandatory. Failing to perform them entails severe risks as an unnecessary asynchronous ventricular stimulation and a pacemaker-mediated tachycardia

**WARNING!**

These tests are mandatory. Failing to perform them entails severe risks as loss of capture, retrograde conduction, induced tachycardia, loss of hemodynamic support

**WARNING!**

If retrograde P-waves are detected out of the rate dependent automatic post ventricular atrial refractory period (PVARP), raise the PVARP manually until such waves are brought back into the PVARP. Failing to perform this procedure may produce a pacemaker-mediated tachycardia

**WARNING!**

A sudden interruption of stimulation may induce asystole before the intrinsic rhythm is restored. Before halting stimulation, set the device on a demand-mode, then gradually reduce the rate below the patient's spontaneous one.

**WARNING!**

Use Pause Button with care because as long as it is pressed the patient has no stimulation (in any case, the button is automatically disabled after 10 seconds).

**WARNING!**

Setting a too long AVD can let the ventricle be stimulated during the vulnerable repolarization period, then triggering ventricular arrhythmia in unstable patients

**WARNING!**

Use the EMERGENCY button only when it is necessary to deliver a high intensity signal (V00). See Sect. 3.4.

**WARNING!**

In the DDI and DDD stimulation modes the signal delivered to the atrium could be detected by the ventricular sensing amplifier. To avoid such a situation it could be necessary to reduce the amplitude of the atrial stimulation and the ventricular sensitivity (raising the corresponding sensing threshold) and/or relocate the electrodes, since they could be too close to each other.

## 5.1 Risks related to electric and electromagnetic disturbances

- An endocavitary or epimyocardial electrode connected to any medical device provides a direct, low resistance current path towards the heart. The risk of inducing menacing ventricular arrhythmia from leaks of alternating currents increases noticeably when a device powered from the mains is connected to the derivation system. Great attention must be paid to ground all devices connected to the mains used in proximity to the patients.
- It is necessary to avoid touching with bare hands metallic parts of the device conductors and letting them come into contact with wet or electrically conductive surfaces.
- A good safety norm is equalizing the potential between the user and the patient. Therefore, it is recommendable to touch the patient in points far from the stimulation derivations.
- Keep all possible source of static electricity far from the stimulation system.
- The radiofrequency generated by an electric scalpel can reach the electrodes in contact with the heart wall, possibly producing arrhythmia or burns on the myocardium. Therefore, such devices have to be used at a sufficiently large distance from the pacemaker and the connected system of electro-catheters.
- Devices generating ionizing radiations (cobalt machines, linear accelerators, etc.) may damage the electronic components in the stimulator, particularly when they accumulate. In case the stimulator is used in such conditions, it is necessary to safely shield the machine and to continuously monitor the patient.
- Not NMR compatible - The temporary pacemaker is not NMR-compatible. Do not take the device to the MR-room.
- After the use of a high frequency medical equipment or a defibrillator in the presence of the device, the latter has to be inspected to be sure that it is still working properly.

## 5.2 Checking the electrical connection

- Make electrical connections following the Instructions for Use at Section 13 "Connection to patient and conditions to use".
- Before any use during heart surgery, check that the clamping jaws on the two bushings of the cable for cardiac surgery (e.g. mod. FIAB F7818/30) close uniformly and firmly on the terminals; in case of doubt do not use them and contact FIAB.

## 5.3 Checking integrity and compatibility

- Before use check that the device and the accessories are intact. Inspect cables and connectors with particular care. In case of doubt, do not use the device and contact FIAB.
- After an accidental fall of the device check that it works properly and in case of doubts contact FIAB before use. If any mechanical damage is evident do not use the device and contact FIAB.
- Before use check that the electro-catheter connector is compatible with the device.

## 5.4 Environment and storage

- For the device maintenance, refer to the IFU at Sect. 11 "Maintenance".
- The environment in which the device is used must comply with the safety norms.
- After storage in an environment not compliant with storage standards concerning temperature and humidity, wait before use until the device stabilizes to adequate temperature and humidity as specified in Table 1 "Technical characteristics", at the end of this document.
- The device must not be immersed in water or any other liquid. Do not pour liquids on the device, though it has been designed to prevent infiltrations.

- Do not place the device where it could be operated, even accidentally, by patients or visitors. Any alteration of the parameters set in the device could have immediate and serious consequences on the patient. Such risks should be reduced as much as possible.

## 6 PERIODIC CONTROLS

The user should control periodically that the devices, its accessories and electric connections are in perfect conditions.

## 7 CLEANING

The device has to be kept clean for hygienic and safety reasons.

## 8 CONNECTIONS TO PATIENT

The device connections can be provided in one of the following configurations: protected bushings, bushings to tighten, Hypertronics connector.

### Monopolar terminals


| Channel             | Symbol            |                   | Terminal labelling  |
|---------------------|-------------------|-------------------|---------------------|
|                     | Positive terminal | Negative terminal |                     |
| Atrial channel      | A+                | A–                | Atrial channel      |
| Ventricular channel | V+                | V–                | Ventricular channel |

### 8.1 Protected bushings

These output bushings can accept 2mm jacks, protected or not.

The protected jack of the electro-catheter or of the extension cable, if any, must possess the following characteristics:

- insulator with internal diameter of 4.5mm and outer diameter of at most 6.2mm (suggested thickness 0.8mm)
- internal contact diameter 2mm

**Figure 1: Protected bushing**

## 8.2 Bushings to tighten

If the device has to be connected to an epicardial needle it is necessary to use either bushings to tighten or a sterile bipolar adapter cable for heart surgery with bushings to tighten. These bushings can receive sockets with diameter between 0.5mm and 2.0mm. For an ideal coupling, it is necessary that the bushel jaws maximal opening is only slightly larger than the diameter of the socket to be inserted.


**Figure 2: Bushing to tighten**

## 8.3 Bipolar terminals

The atrial and ventricular bipolar pairs endings are separated and indicated with the respective letters.

### Hypertronics Connector

This type of connector is frequently used for external pacemakers

**Figure 3: Connection diagram**

# 9 TECHNICAL CHARACTERISTICS

|                                                                    |                                                                                                                                                                                                                                                                                                                                                                                                                                         |
|--------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Producer</b>                                                    | FIAB SpA Via Costoli, 4 - 50039 Vicchio - Firenze (ITALIA)                                                                                                                                                                                                                                                                                                                                                                              |
| <b>Model</b>                                                       | 1797                                                                                                                                                                                                                                                                                                                                                                                                                                    |
| <b>Class and Type</b>                                              | Portable device, internally powered                                                                                                                                                                                                                                                                                                                                                                                                     |
| <b>Applied part<br/>Degree of<br/>liquid/solid<br/>penetration</b> | CF type, protected from defibrillator discharge<br>Complying to particular standard EN 60601-2-31:2020                                                                                                                                                                                                                                                                                                                                  |
| <b>Conditions of use</b>                                           | Device for continuous operation                                                                                                                                                                                                                                                                                                                                                                                                         |
| <b>Compliance to<br/>safety standards</b>                          | EN 60601-1:2007/ A1:2014 Electric MDs - Part 1: General safety standards<br>EN 60601-2-31:2020 Electric MDs - Part 2: Particular safety standards for internally powered external cardiac stimulators<br>EN 60601-1-2:2015 "Medical electrical equipment – Part 1-2: General requirements for basic safety and essential performance – Collateral Standard: Electromagnetic disturbances – Requirements and tests (IEC 60601-1-2:2014)" |
| <b>Intended use</b>                                                | endocavitary temporary single and dual chamber stimulation                                                                                                                                                                                                                                                                                                                                                                              |
| <b>Operation<br/>conditions</b>                                    | Temperature: from +10C to +40 C<br>Relative humidity: from 30% to 70% RH not condensed<br>Atmospheric pressure: from 700 mbar to 1060 mbar                                                                                                                                                                                                                                                                                              |

|                                                            |                                                                                                                                                                                                            |
|------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Storing and transportation conditions (4 weeks max)</b> | Temperature: from 0C to +50 C<br>Relative humidity: from 20% to 80% RH not condensed<br>Atmospheric pressure: from 500 mbar to 1060 mbar<br>Packaging suitable to prevent shock induced mechanical damages |
|------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|

| Rev. | Revision Status / Description | Date       |
|------|-------------------------------|------------|
| A    | First Issue                   | 22/11/2023 |
| B    | Second Issue                  | 12/11/2024 |
| C    | Third Issue                   | 15/06/2026 |
| D    | Fourth Issue                  | 10/09/2026 |

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