# 1748

## **SINGLE-CHAMBER EXTERNAL CARDIAC STIMULATOR** **STIMOLATORE CARDIACO ESTERNO MONOCAMERALE**

# INFORMATION FOR SAFETY AND PERFORMANCE

## (EN)

### 1. DEVICE DESCRIPTION

FIAB 1748 external pacemaker is a compact, temporary single-chamber cardiac stimulator to be used in case of pathological alterations of the heart pulse spontaneous generation or conduction. Powered by an internal battery, it is connected to electro-catheters for temporary atrial or ventricular stimulation in clinical environment. Various stimulations modes are available, along with the high-rate atrial BURST and RAMP functions.

The stimulation mode can be selected via a quick menu together with pulse rate and amplitude, and the input signal sensitivity (sensing threshold); the high-rate stimulation modes BURST and RAMP can be selected, as well as the BURST rate and the maximal RAMP rate.

The FIAB 1748 principal characteristics are:

- Stand by, pre-stimulation modality for parameters and modes selection
- synchronous and asynchronous stimulation;
- system driven by two microprocessors, allowing high precision in selecting the stimulation parameters and a complete and continuous control of the device functioning;
- low consumption liquid crystals display;
- EMERGENCY function for a quick launch of asynchronous stimulation with safety parameters;
- BURST and decreasing RAMP modes for high frequency stimulation in the treatment of atrial tachycardia;
- backup batteries allowing main battery replacement without interrupting stimulation;
- easy and safe connection to temporary electro-catheters for transvenous epicardial or myocardial stimulation, either direct or by means of a separate patient cable or with an adaptor, when needed.

#### 1.1. Safety functions

- Self-test function;
- On-line check of batteries charge level, signaling exhaustion both of the main and of the backup battery
- Continuing stimulation during battery replacement;
- Keyboard lock preventing accidental parameters modification;
- Two-step switching off preventing accidental deactivation of the device;
- Anti-runaway protection preventing stimulation rate to go out of control;
- Protection against defibrillator discharge;
- Protection against electrostatic discharge;
- Protection against electromagnetic interferences;
- Continuous check of electrical connections, of load impedance and of stimulating current;
- Potentially dangerous situations trigger alarms to the user by means of acoustic signals and corresponding warnings on the display

#### 1.2. Main operation features

- Single chamber stimulation modes: S00, SST, SSI (S = Single Chamber, Atrium or Ventricle)
  - S00 asynchronous mode;
  - SSI synchronous mode;
  - SST triggered mode;
- Pacing and sensing LEDs, showing that electric signals are being sent (pace) or received (sense) by the device;
- Low battery messages warning that the main battery has to be replaced;
- Low battery messages warning that the back-up battery has to be replaced;
- Display devoted to the setting of stimulation characteristics (amplitude, rate, sensitivity) viewable and easily adjusted with knobs;
- Alarm and warning messages shown in the monitor lower band;
- Keyboard lock/unlock key protecting against unintentional parameter changes;

- Pause key – absence of stimulation – suspending stimulation and allowing inspection of the patient's own intrinsic rhythm;
- Emergency key – starts single chamber emergency asynchronous stimulation (S00).

### **1.3. Indications**

Temporary pacemaker is used in a hospital setting for temporary single chamber cardiac pacing in case of pathological alterations of the heart pulse spontaneous generation or conduction:

- Cardiac arrest or complete heart block
- sinus bradycardia

### **1.4. Intended users**

The 1748 single-chamber stimulator is intended to be used by highly qualified medical personnel: specialists in cardiology or cardiac surgery, experts in cardiac stimulation. Before starting to use the device, train medical personnel on the functionality and use of the temporary pacemaker, based on needs and hospital procedures.

### **1.5. Environment of use**

The DM must be used in hospital environments: cardiology, emergency room and cardiac surgery. The temporary pacemaker must be used in an environment where continuous monitoring of the patient is guaranteed, in order to ensure optimal functioning of the device and guaranteeing that the patient receives the most adequate treatment.

### **1.6. Essential Performance**

The essential performances of the MD can be summarized in two fundamental functions.

- Pacing: the MD must stimulate the heart tissue (with amplitude and frequency that can be set by the user and in the chosen mode - synchronous or asynchronous)
- Sensing: the MD must recognize spontaneous cardiac activity and synchronize with it if the device has been set with ON DEMAND functionality (synchronous).

### **1.7. Clinical Benefit**

Temporary cardiac stimulation is used in the treatment of heart disorders to:

- Heartbeat recovery, keeping the patient alive until spontaneous recovery or until the implantation of a permanent pacemaker
- Maintenance of the correct heart rate in case of pathological alterations of the spontaneous generation of cardiac impulse, or stabilization of the patient from a hemodynamic point of view, i.e. guaranteeing the correct flow of blood into the patient

### **1.8. Intended purpose**

A temporary pacemaker is used in connection to a system of electro-catheters, so to provide a temporary single- or dual-chamber stimulation in a hospital environment by trained personnel. The temporary pacemaker can provide either synchronous (on demand) or asynchronous, short-term stimulation for therapeutic, prophylactic or diagnostic purposes, in particular:

- support, management, assessment of a patient with cardiac complications during invasive procedures
- patient support after heart surgery

### **1.9. Intended patient population**

Patients who need short-term, synchronous or asynchronous cardiac electro-stimulation for therapeutic, prophylactic or diagnostic purposes.

## **2. BATTERY**

The device is powered by an internal battery. Only 9 V batteries of type 6LR61 (identification in accordance

with international standard IEC 60086) may be used. Batteries of different physical sizes or rechargeable batteries can cause the temporary pacemaker to function erratically, fail to pace the patient, or damage the temporary pacemaker. The service life of rechargeable batteries is unpredictable and their use may result in sudden termination of stimulation. Before each use, especially if prolonged, it is strongly recommended to replace the battery. This replacement then becomes indispensable if the battery change alarm appears on the stimulator. Only use new, charged batteries that have not passed their expiration date. If the stimulator is stored for an extended period, the battery must be removed from its compartment.

### 3. POTENTIAL COMPLICATIONS – RESIDUAL RISKS

Concerning stimulation in general:

Temporary external stimulation may cause complications, among which:

- asystole after sudden pacing interruption (e.g. caused by unadverted accidental disconnection of the patient cable, of the electrocatheter or by wrong settings);
- inhibition or decay because of intense electromagnetic interferences;
- onset of tachyarrhythmia or worsening of an already present one;
- dependence on the pacemaker.

Concerning the specific device:

- muscular contractions, possibly painful;
- stimulation irregular or absent due to electro-catheter breakage or dislocation;
- loss of the signal capture due to a wrong regulation of the stimulation parameters (ineffective stimulation);
- tachycardia induced by a wrong regulation of the stimulation rate (inappropriate stimulation);
- loss of synchronism or stimulation inhibition due to a wrong regulation of the sensing parameters (Undersensing or Oversensing).

### 4. CONTRAINDICATIONS

No contraindications are known advising against the use of temporary stimulation to control heart functions. However, patient's age and clinical conditions may limit the choice of stimulation modes and parameters and of the most adequate catheter system.

#### Atrial stimulation

Atrial stimulation is contraindicated in the presence of atrial fibrillation or flutter. The single chamber atrial stimulation is contraindicated in the presence of imperfect A-V conduction.

#### Asynchronous stimulation

Asynchronous stimulation is contraindicated in the presence of intrinsic heart rhythm.

#### High-rate atrial stimulation therapy

High-rate stimulation therapy can only be used in atrium. If used in ventricle can have lethal consequences.

#### Competitive stimulation

Temporary stimulation is contraindicated in the presence of another stimulating system. Do not use the temporary pacemaker while another stimulation system is working. This may result in a competition between the two systems. In that case the temporary pacemaker can be inhibited or the stimulation can take place in an asynchronous way. As a consequence, a stimulation during a T-wave or a pacemaker mediated tachycardia could intervene.

### 5. WARNING AND PRECAUTIONS

The necessary precautions for the correct use of the single-chamber stimulator are given below. Other indications are described in the IFU.

The manufacturer is held responsible for the effects of the safety, reliability and performance of the appliance only if these warnings are respected.

- **Who can use it and when.** The 1748 stimulator is intended to be used by highly qualified medical

personnel: specialists in cardiology or cardiac surgery, experts in cardiac stimulation. The DM must be used in a hospital setting: cardiology, first aid and heart surgery.

- **Read instructions.** The device must be used in accordance with the operating instructions. It is necessary to read the IFU entirely, if possible in front of the appliance which is switched on and running, in order to fully understand its functions.
- **ECG Monitoring.** It is essential to check the effectiveness of the stimulation and the correct setting of the stimulation parameters using an ECG monitor. During the use of the stimulator it is essential to monitor the patient's ECG. In case of prolonged use of the stimulator, periodic checks by the doctor of the sensing and pacing thresholds are necessary. For emergency situations, you should also have a defibrillator ready to use during the following procedures: evaluating the pacing and sensing thresholds, connecting the pacemaker to the pacing leads and making adjustments, as well as during high-rate Burst or Ramp atrial pacing therapy.
- **Sensing threshold.** An appropriate setting of the sensing threshold is essential for the mode SSI to work correctly. If it is too high the device may ignore the patient's intrinsic cardiac activity, triggering unnecessary stimulation (this is commonly called 'Undersensing'). A too low threshold is not recommended either, since irrelevant activities (like e.g. muscle tremors) or slight disturbances could be erroneously interpreted as spontaneous cardiac activity, inducing continuous inhibitions of synchronous stimulation; moreover, overlapped disturbances could hide the physiological signal making it interpreted as a noise with the consequence that the device will ignore it (these are the two aspects of the so-called 'Oversensing').
- **Asynchronous stimulation.** The asynchronous stimulation mode (S00) may interfere with the patient's spontaneous activity and therefore is normally not recommended. It is recommended only in the absence of any spontaneous activity.
- **Accessories.** For the parts in contact with the patient follow the IFUs provided with them (see IFU sect. 14).
- **Sensing.** Since the value of sensitivity determines the minimal signal detectable by the device, it is necessary to set it to at least half of the determined sensing threshold (see the IFU sect. 6 Determination of the sensing threshold). This guarantees a sufficient safety margin for the sensing activity. However, pay attention not to select a too low value otherwise far field signals can be sensed (like the R and T-wave when sensing in atrium or the P-wave when sensing in the ventricle), causing the wrong inhibition of the stimulation pulses.
- **Safety Margins.** Take adequate safety margins both for sensing and for pacing, otherwise the stimulation provided by the device could be inappropriate.
- **Sensing threshold of implanted devices.** Do not use the temporary stimulator to determine the sensing threshold of implanted systems. During the implantation of a permanent pacemaker, FIAB recommends the use of a threshold analyzer (PSA e.g.: FIAB 2307).
- **High-rate atrial Burst or RAMP stimulation therapy.** The use of BURST and RAMP high-rate stimulation is basically indicated for the treatment of atrial synchronized tachyarrhythmia. The Burst therapy can only be used in atrium, its use in the ventricle may induce lethal arrhythmia. The same is true for the Ramp therapy.
- **High-rate atrial Burst or RAMP stimulation therapy.** Atrial high-rate stimulation may induce fibrillation, tachycardia, or accelerate an already existing tachycardia. Atrial high-rate stimulation must be performed in a constantly controlled clinic environment and it has to be monitored, checking the ECG and the blood pressure. A defibrillator must be readily available, possibly already connected to the patient with adhesive pads.
- **Sensitivity threshold test.** Perform the sensitivity threshold test to determine the most adequate sensing threshold. Among the clinical risks associated to skipping this operation are possible asynchronous ventricular stimulation and a pacemaker mediated tachycardia
- **Pacing threshold test.** Perform the stimulation threshold test to determine the most adequate stimulation amplitude. Among the clinical risks associated to skipping this operation are loss of capture, retrograde conduction, induced tachycardia and loss of hemodynamic support.
- **Stimulation interruption.** Sudden stimulation interruption may produce asystole. Before stopping a stimulation set the pacemaker on SSI mode, then gradually reduce the stimulation rate taking it below the patient intrinsic rate.
- **Pause button.** Use the Pause button with care, since, as long as it is pressed, the patient is not stimulated (the system disables it anyway after 10 sec.).
- **Emergency button.** Use the Emergency button only when it is necessary to supply a high output asynchronous stimulation (S00).

### **Risks connected to electric or electromagnetic disturbances**

In general, an endocavitary or epimyocardial electrode connected to any electromedical device constitutes a direct current path, with low resistance, towards the myocardium. The risk of inducing menacing ventricular arrhythmias resulting from AC leakage greatly increases when line-powered equipment is connected to the shunt system. Great care should be taken to ground all mains powered equipment used in the vicinity of the patient.

- Do not touch the metallic part of the terminals with bare hands, nor put them in contact with wet or electrically conduction surfaces.
- Since it is always good practice to equalize the potential between the User and the Patient, the User is recommended to touch the Patient at points away from the stimulation leads.
- Keep all possible sources of static electricity away from the system.
- The presence of a lead in contact with the heart wall can allow the radiofrequency generated by an electrosurgical unit to reach the heart wall itself, producing possible arrhythmia or myocardial burns. Consequently, these units must be used remotely from the pacemaker/lead system.
- Equipments generating ionizing radiation (cobalt machines, linear accelerators, etc.) can damage the electronic components contained in the stimulator, especially following the accumulation of large quantities of radiation. In case of use in the above conditions, effective shielding of the stimulator and continuous monitoring of the patient are required.
- **Not MRI compatible.** The temporary pacemaker is not MRI compatible. Do not take it to the MRI room.
- The device must be inspected after it has been used with high frequency surgical devices or defibrillators.
- The use of accessories and cables other than those specified or supplied by the manufacturer of this equipment (see section 14 of IFU) could result in increased electromagnetic emissions or decreased electromagnetic immunity and could result in improper operation.
- Portable radiofrequency communication equipment (including peripherals such as antenna cables and external antennas) must be used at a distance not less than 30 cm (12 in) from any part of the 1748 device, including the cables specified by the producer. Otherwise, the device could behave improperly.

### **Checking electrical connections**

- Make connections strictly following the instruction in IFU at sect. 13 "Connection to patient and mode of use". Warning: for more safety it is advisable to connect the patient terminals to the bushels of the extension cable before connecting the corresponding jacks to the stimulator.
- Before each new use in cardiac surgery, check that the clamping jaws of the two cable bushings for cardiac surgery (e.g.: model FIAB F7818/30) close uniformly and firmly on the terminals; if in doubt, do not use them and contact **FIAB**.

### **Checking integrity and compatibility**

- Before each use of the device and of its accessories check their integrity. Particular attention must be devoted to inspecting cables and connectors. If in doubt, do not use them and contact **FIAB**.
- After a fall, check that the device works correctly and if in doubt contact **FIAB** before use. In the presence of mechanical damages do not use the device and contact **FIAB**.
- Before each use, check the compatibility of the lead connector with the connection system to the device.

### **Maintenance warnings**

- **Opening the device, accessing its internal parts and repairs must be carried out by FIAB authorized personnel.** Do not modify the temporary pacemaker. This could compromise its effectiveness with negative consequences for the patient safety.

### **Use and storing environment**

- For the device maintenance refer to IFU at sect. 11 "Maintenance".
- The environment in which the device is employed must comply with the safety standards.
- If the device has been stored in a not appropriate environmental temperature and humidity, wait before using it, allowing it to stabilize at adequate temperature and humidity conditions, as indicated in the "technical characteristics" table at the end of this document.
- The device should not be submerged in water of other liquids. Do not pour liquids on the temporary pacemaker. The device has been designed to reduce infiltrations to the minimum, but they could take place anyway.
- To not put the temporary pacemaker where it could be operated, even accidentally, by the patient. The alteration of the set parameters could have direct and serious effect on the patient. The pacemaker must be put in an area minimizing the risk of tampering by non-authorized persons (patients or visitors).

- The temporary pacemaker is not suitable for use in the presence of flammable anesthetics mixed with air, oxygen, or nitrous oxide (N<sub>2</sub>O).

**Do not sterilize**

- The product is not sterile and cannot be sterilized: if the device is put in the sterile area, hygienic conditions are not guaranteed.
- To avoid damages, the device must not be sterilized. To clean the device and its accessory follow the indications in the section 11.2 of the IFU.
- Accessories provided with the device are not sterile.
- Pay attention not to contaminate the patient cable sockets with blood or other body fluids.

## 6. CHECKS

The user must periodically check the functionality of the device, of its accessories and on the connections state. If any malfunction or loosening of the connections is detected, the FIAB Technical Assistance Service must be contacted.

Moreover, the device must receive an ordinary maintenance intervention at least once a year at FIAB or by FIAB authorized personnel. Ordinary maintenance includes the following operations:

- replacement of back-up battery;
- functional verification;
- electric safety verification;
- final test.

## 7. CLEANING

It is essential that the user takes care of cleaning operations requested by ordinary hygienic needs. Clean the device removing all visible traces of impurities or blood.

## 8. CONNECTION TO PATIENT

Concerning the input/output terminals, the device can be provided in one of the following configurations: protected sockets, bushings to tighten, Hypertronics connection.

### 8.1. Monopolar Terminals


Terminals polarity is specified on the device side hosting connections.  
Terminals (+) are red; Terminals (-) are black.

#### Protected sockets

The outlet protected sockets can accept 2mm protected or unprotected pins. To prevent connection problems, be sure that the electrocatheter plugs are compatible with the outlet sockets of the device, so that both electrical and mechanical stability are guaranteed.


The electrocatheter protected pin (or the one of the extension cable) must have the following features:

- insulator with 4.5mm minimal internal diameter and 6.2mm maximal external diameter (recommended thickness 0,8 mm).
- 2mm internal contact diameter.

Respect connection polarities.


Insert the touch-proof pins fully, beyond the first resistance encountered, in order to obtain complete coupling between the pin and the metal insert of the bushing.


#### Bushings to tighten

When connecting to an epicardial needle it is necessary to use the bushings to tighten version or a sterile bipolar adapter cable for cardio surgery (e.g. F7818/30) equipped with bushings to tighten. The bushing may accommodate terminal with diameter from 0.5mm to 2.0mm. For an ideal coupling it is necessary that the opening of the jaws of the bush has a diameter equal to or slightly greater than that of the terminal to be inserted. The diameter of the opening increases (up to 2 mm) by rotating the bushing protection counterclockwise, and decreases (up to 0.5 mm) by rotating it in the opposite direction.


The epimyocardial needle is commonly characterized by a diameter of about 1mm. Before inserting it into the bushing, it is advisable to visually check that its opening is adjusted in the appropriate way. Furthermore, in order not to compromise the efficiency of the clamping system, it is advisable to insert the needle keeping it in line with the axis of the bushing and then to tighten it well by turning the bushing clockwise, to guarantee effective electrical and mechanically safe contact. Finally, it is advisable to protect (for example with FIAB 'Myoprotection') the residual metal part of the needle that may protrude from the insulating coating of the sleeve.


## 8.2. Bipolar Terminals

### Hypertronics Connector

The Hypertronics connector is widely used for external pacemakers (particularly on pacemaker by Medtronic).

Check compatibility between the catheter and the device connectors.  
 A tab forces the insertion of the connector in a predefined position, preventing accidental insertion in the inverted configuration. Fully insert the connector until you hear the click of the tab which locks the connector in its housing.

## 9. TECHNICAL CHARACTERISTICS

Producer	FIAB SpA Via Costoli, 4- 50039 Vicchio - Firenze (ITALIA)
Model	1748
Class and Type	Portable device, internally powered
Applied part	CF type, protected from defibrillator discharge
Employment conditions	Device for continuous operation
Compliance to safety standards	EN 60601-1:2006+A2:2021 Medical electrical equipment - Part 1: General requirements for basic safety and essential performance EN IEC 60601-2-31:2020 Medical electrical equipment - Part 2-31: Particular requirements for the basic safety and essential performance of external cardiac pacemakers with internal power source EN 60601-1-2:2015+A1:2021 Part 1-2: General requirements for basic safety and essential performance - Collateral Standard: Electromagnetic disturbances - Requirements and tests MDD 93/42/CEE
Intended use	Endocavitary temporary single chamber stimulation
Operation conditions	Temperature: from +10°C to +40 °C Relative humidity: from 30% to 70% RH not condensed Atmospheric pressure: from 700 mbar to 1060 mbar
Storing and transportation conditions (4 weeks max)	Temperature: from 0°C to +50 °C Relative humidity: from 20% to 80% RH not condensed Atmospheric pressure: from 500 mbar to 1060 mbar Packaging suitable to prevent shock induced mechanical damages

Rev.	Revision Status / Description	Date
A	First Issued	12/09/2023
B	Second Issued	12/11/2024

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