SAFETY AND PERFORMANCE SUMMARY

ADAPTERS FOR ELECTROSURGERY PENCILS

EN

1. PRODUCT DESCRIPTION

Pencil adapters are non-sterile reusable devices, class I according to EU Regulation 2017/745 (MDR), intended for adapting the pencils to the electrosurgical unit connector.

2. INDICATIONS AND USERS

The adapters are devices of a very limited length (max 15 cm) and built in a single body. They are equipped on one side with a connector that fits into the electrosurgical unit, specific according to the manufacturer and model of the latter. On the other side there is a female or male connector which allows you to adapt the connection of the electrosurgery pencil or the electrode you intend to use.

The device is intended for use by individuals who have had training in the use of electrosurgical units.

A preliminary procedure is recommended before using the device. During the above procedure always select the lowest possible energy level. If the clotting capacity of the electrode is less than normal, do not increase the high frequency output without having previously carried out the following checks:

- the correct positioning of the neutral plate;
- the correct insertion of the cables and their connectors;
- the correct activation of the operating keys (manual or pedal);
- there is no insulation damage of the cables and connectors.

3. INTENDED PURPOSE

Adapters are reusable devices designed to adapt the different types of electrosurgical units to the electrosurgery pencils.

4. ESSENTIAL PERFORMANCE

The adapters allow the transmission of high frequency energy from electrosurgical unit to the patient.

5. CONTRAINDICATIONS

The products should not be used if:

- the patient has a pacemaker unless prior cardiological advice;
- there is visible damage to the adapter, its insulation, the body and/or the connection cable;
- the preliminary procedure (see ch. 2) gave a negative result;
- the electrode is not inserted firmly into the distal part of the product.

6. POSSIBLE COMPLICATIONS

Potential complications can be identified as burns caused by a high current density in the patient's tissue or by heating flammable fluids or gases; the causes may be the following:

the patient was inadvertently placed in contact with electrically conductive parts;

there was direct contact between the cables and the patient's skin resulting in a capacitive effect;

- combustion of flammable disinfectant agents;
- combustion of flammable narcotic gases.

6.1. Residual risk

There are no residual risks if the use of the device is in accordance with the instructions for use attached to it.

7. PRECAUTIONS AND WARNINGS

These precautions should be followed to reduce the complications listed above.

7.1. Sterility

MULTI-USE DEVICE – NON STERILE

The devices MODELS F4830; F4834; F4835; F4836; F4837; F4838; F4831; F4832 and F4833 can be sterilized with ethylene oxide (EtO)

7.2. Equipment and accessories

The adapters must be used with electrosurgical units complying with EN 60601-1, EN 60601-1-2 and EN 60601-2-2 standards. Before use, check the instructions for use for any restrictions.

7.3. Electrical and electromagnetic safety

Pacemakers can be damaged by an electrosurgical current. Do not submit a pacemaker patient to electrosurgical current without first consulting a cardiologist.

7.4. Performance and use

The procedures for fitting and removing the adapters, as well as cleaning, must only be carried out with the generator switched off.

Always carry out a functional check before use.

Reject defective products.

Make sure that the connection of the pencil plug to the generator is stable and that the metal part of the connector is fully inserted.

Keep the connected adapters and pencils away from flammable materials: their accidental activation and the heat could cause fires.

When the products are not used, store them in suitable containers.

HF electrosurgical generators, paddles and accessories to be used with these products must comply with applicable regulations.

For the use of the generator, pencils and neutral plates, refer to the instructions for use provided by the manufacturer.

The products must be used by qualified healthcare personnel.

Check the compatible adapter with the electrosurgical unit and the pencil used.

Do not force liquids or air into the adapter.

Do not use ultrasonic baths.

The use of mechanical washing systems can reduce the life of the product.

If the device is sterilized, the validation of the process is under the responsibility of the user.

For ethylene oxide sterilization, clean products must be packaged in a suitable gas-permeable bag. It is recommended to respect an adequate aeration time to allow total dispersion of EtO residues.

The product must be stored in its original packaging in rooms characterized by environmental conditions, temperature and relative humidity, specified on the label on package.

Overlapping weights on packages could damage the product.

Rev.	Revision Status / Description	Date
A	First Issued	23/02/24

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