Connection and extension cables

INFORMATION FOR SAFETY AND PERFORMANCE

(EN)

1 INTENDED USE

1.1 Single use (sterile) and reusable extension cables

Devices intended for connection between a temporary epimyocardial electrode or bipolar Endocavitary electrode and an external cardiac stimulator.

1.2 Connection cables

Connection cables may be used as a connection between the defibrillator and the electro catheters during endocavitary defibrillation procedures.

2 PRECAUTIONS AND WARNINGS

2.1 Single use (sterile) and reusable extension cables

- **Check the package before using the product. If wet or torn its sterility is not guaranteed (*only for single use-sterile*).**
- **Sterilize before use (*only for Reusable*).**
- The devices that are connected to the cables must comply with the law in force.
- Only use devices that are suitable for direct cardiac application (like CF).
- NB: prior to connection, read the instructions manual of the equipment to be used.
- All electric equipment supporting the patient must be grounded via equipotential node.
- The devices must be connected and operated by qualified staff.
- Prior to use, check the integrity of the cable coating, of the plastic parts and of the connections; do not use the product if it is damaged or faulty.
- The cable must be connected to the equipment as the first connection of the system. The extension cable must be connected to the electrode as the last connection of the system.
- Before each procedure, carefully check that the plugs have been properly inserted and that all connections are solid, pulling them gently to check eventual movements.
- Work in a sterile environment.
- Make sure that a defibrillator is ready for use.
- The product life is reported on the label.

WARNING!

- BIO" suffix models can be used only with PACEMAKER BIOTRONIK Mod. Reocor S.
- "MED" suffix models can be used only with MEDTRONIC PACEMAKER Mod. 5392 and Mod. 5388 (Dual Chamber) or 5348 (Single chamber).
- "30" suffix models can be used only with devices fitted with sockets complying with EN60601-1 and E-DIN 42802-2 regulations.
- When using the product, follow all the precautions and warnings reported below, especially those concerning models with "33" suffix since the connections that are far from the patient are uncovered once the protection caps are removed.
- "OSY" suffix models can be used only with OSYPKA devices with screw sockets.

2.2 Connection cables

- The product supplied is not sterile
- Devices connected to cables must all be complying with the regulations in force and of CF type.

- NB: before making connections consult the instruction manual pertaining to the instrument to be used.
- All the electrical instruments must be opportunely grounded
- Medical devices must be connected and activated by qualified personnel only.
- Prior to use check the integrity of the cable sheath, the plastic parts and the connections. Do not use the product if there are any signs of defects or damage.

Only for F5414 series

- Fully insert the male connectors into the female connectors in order to ensure that no sneak electric currents occur during the delivery of the discharge from the defibrillator. The same precaution applies for all the connections of the F5414/XX cables to the devices.
- Before each operation carefully check the correct insertion and hold of all the connections.

Only for F7840 series

- Do not wash the leads with a solvent; if necessary, during the procedure, use a sterile physiologic solution. Dry well. Even during cleaning and drying it is necessary to stay within the temperature limits indicated on the label (0°C-50°C).
- Before the procedure check the effective insertion of every pin and connection by gently tugging them, to feel any movement.
- Whenever using an ECG monitor and/or a radiofrequency generator for positioning the electrodes, it is necessary for the equipment and the cable to be insulated. Loss of current by the monitor must not exceed 10 µA for intracardiac electrodes.

3 POTENTIAL COMPLICATIONS

3.1 Single use (sterile) and reusable extension cables

An incomplete blockage of the electrodes to the cable can cause an intermittent or ineffective stimulation and can provide incomplete protection against dispersion currents, with the ensuing risk of dangerous ventricular arrhythmias.

3.2 Connection cables

An imperfect socket-plug connection could represent an electrocution hazard for the operator if he/she accidentally makes hand contact during the delivery of the defibrillator discharge.

Rev.	Revision Status / Description	Date
A	First Issued	13/01/2025

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