EURODEFIPADS®

DISPOSABLE MULTIFUNCTION ELECTRODES

INFORMATION FOR SAFETY AND PERFORMANCE

(EN)

INTENDED USE

EURODEFIPADS disposable multifunction electrodes are indicated for:

- transthoracic external defibrillation,
- transthoracic synchronized cardioversion,
- transthoracic ECG Monitoring,
- temporary transthoracic cardiac pacing (non-invasive).

The device is intended to be used in non-sterile environment by qualified healthcare personnel and/or, if applicable, by people trained for CPR (Cardio-Pulmonary Resuscitation) and for the use of AED (Automatic External Defibrillator).

Adult models are addressed to patients that weigh more than 25 kg.

Paediatric models, REF with P suffix (F7xxxP), are addressed to children that weigh less than 25 kg.

Adult/paediatric models, REF with U suffix (F7xxxU), are addressed to patients that weigh more than 10 kg.

CONTRAINDICATIONS

- The use of FIAB adult disposable electrodes is generally contraindicated in patients younger than 8 years old (weight less than 25kg), but can be used if the size of the chest allow it, taking care that one pad do not come in contact with the other one and following the operating instructions of the defibrillator for energy to be delivered.
- The use of FIAB disposable multifunction electrodes in versions for adult or adult/paediatric is contraindicated in children aged less than 12 months (weight less than 10kg). In this case use a pediatric model.
- The use of FIAB pediatric disposable electrodes is contraindicated in patients older than 8 years old (weight more than 25kg). In this case use a model for adult or adult/pediatric.
- Do not apply on skin that shows signs of irritation or injury.

PRECAUTIONS AND WARNINGS

- Use the device only with compatible defibrillators. Check the compatibility on the package.
- Read the instructions for use of the defibrillator, with particular attention on the placement of multifunction electrodes, their polarity and the power to be supplied.
- In paediatrics and for some models of automatic defibrillators the use of specific reducing power devices or the adoption of special precautions may be required. Always pay a special attention to energy levels set on the defibrillator and that can be delivered to the paediatric patient.

ATTENTION

Do not use paediatric multifunction electrodes marked with this symbol with automatic defibrillators.

Paediatric multifunction electrodes marked with the symbol shown beside are indicated for use with automatic defibrillators.

- The electrode choice should be based on the evaluation of chest size and weight of the patient. Paediatric electrodes used beyond the specified energy limit may cause also major skin burns; on the contrary the extended active surface of adult electrodes may jeopardize the therapy when used for pediatric treatment.
- After an extended period of transthoracic pacing the ability to detect the evoked ECG signal can be reduced. In this case it is necessary to provide by a separate set of ECG electrodes for the collection of the evoked signal.
- Replace the multifunction electrodes after 24 hours from their application on the patient's skin.
- Do not use the electrodes if they exceed the expiration date reported on the package.
- The traceability data and the device expiration date are detailed only on the package: save the package or write REF and LOT # for possible references to the used electrodes.
- Check that the packaging is intact: do not use the product otherwise.
- Open the package of the multifunction electrodes only before they use. The adhesive pads must be attached on the skin of patient immediately after the removal of the protective covering.
- Do not use the multifunction electrodes if the gel is removed from the support or if it is ripped, torn or dry. Any discoloration localized on gel or on conductive foil does not affect the functionality of the product.
- Do not use the multifunction electrodes if during removal of the protective coating the product is damaged (eg. the insulating coating of the contact has detached or there are tears in the foam support and/or in the electrode).
- Do not bend, cut nor squash the adhesive pads.
- Do not use the multifunction electrodes if the connector, the cable or the clips appear to be damaged.
- Check on the operating instructions of the defibrillator at which safety distances the devices (surgeon's electric knife, RF ablaters, diathermy equipment, mobile phones, etc.) that emit strong electromagnetic interferences must be placed. Position the defibrillator/electrode system at a distance of at least one and a half times the recommended separation distances.
- The electrodes and their cable are made of ferromagnetic materials, therefore they must not be used in presence of the intense magnetic field supplied from a magnetic resonance imaging (MRI).
- To prevent accidental damage from electric shock, ensure that during discharge operators are not in contact with the electrode pads, with the patient, or with conductive parts close to the patient.
- When defibrillators are used near oxygen sources or other flammable gases, use extreme care to avoid risk of fire or explosion.
- The device is not sterile. Neither disinfect nor sterilize it.
- The device is disposable - do not reuse. The reuse can imply: material alterations and the loss of the device functional characteristics.

POTENTIAL COMPLICATIONS

There are no complications related to the use of multifunction electrodes.

WARNING. Defibrillator discharge may cause irregularities in the operation of an implanted pacemaker/defibrillator; apply the multifunction electrodes at a distance of at least 8cm from the said devices. After defibrillator discharge check its operation.

WARNING. If the chosen energy level is insufficient the success of therapy may be jeopardized. On the contrary, levels higher than necessary may modify the enzyme structure without actual evidence of myocardial damage.

Rev.	Revision Status / Description	Date
A	First Issued	26/09/2022
B	Translation order reversal	13/10/2022
C	Introduction of form 99501185	03/07/2023
D	Intended patient population specified in more detail for each type of electrode	18/02/2026

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