
4429

SINGLE-CHAMBER EXTERNAL CARDIAC STIMULATOR **STIMOLATORE CARDIACO ESTERNO MONOCAMERALE**

INFORMATION FOR SAFETY AND PERFORMANCE

(EN)

1. DEVICE DESCRIPTION

External 4429 Easypace pacemaker is a compact, single-chamber, temporary heart stimulator to be used in the pathologies of spontaneous generation and conduction of impulses.

The unit's main characteristics are:

- synchronous and asynchronous modes of operation
- fully microprocessor operated control allowing a high precision setting of stimulation parameters and a complete and continuous check for correct operation of the unit
- back-lighted LCD display
- AUTO and EMERGENCY keys for a quick start-up of both synchronous and asynchronous stimulation with safety parameters
- stimulation capability with quadruple rate and decreasing ramp to reduce tachycardia events
- stand-by batteries for replacement of main batteries without interrupting the stimulation
- a large number of safety functions
- easy and safety connection of Patient cable or lead to sockets
- protected sockets coloured black for the – and red for the +; these are compatible with protected pins fitted with Ø 2 mm electrical- contacts

1.1. Safety functions

- self-test at power on;
- on-line control of the batteries level;
- turn-on inhibition with low battery charge;
- check for run-away and lead connection;
- protection against defibrillation discharges and electromagnetic interference
- Dangerous situations are signalled to the operator by an acoustic alarm and a suitable display message

1.2. Main operation features

The operating characteristics of the temporary pacemaker are as follows:

- single chamber stimulation mode
 - V00 asynchronous mode;
 - VVI synchronous mode;
- Pacing and sensing LEDs, showing that electric signals are being sent (pace) or received (sense) by the device;
- Low battery messages warning that the main battery has to be replaced;
- Low battery messages warning that the back-up battery has to be replaced;
- Display devoted to the setting of stimulation characteristics (amplitude, rate, sensitivity);
- Alarm and warning messages shown in the display;
- Pause key – absence of stimulation – suspending stimulation;
- Emergency key – starts single chamber emergency asynchronous stimulation (S00).

1.3. Indications

Heart stimulator 4429 Easypace is a versatile unit easy to use for single-chamber temporary intracardiac stimulation. In particular it is recommended in:

- extreme bradycardia;
- paroxysmal atrial-ventricular block;
- interruption of paroxysmal tachycardia from atrial flutter;
- atrial-ventricular blocks of various types prior to implantation or replacement of a definitive pacemaker

1.4. Intended users

The 4429 single-chamber stimulator is intended to be used by highly qualified medical personnel: specialists in cardiology or cardiac surgery, experts in cardiac stimulation. Before starting to use the device, train medical personnel on the functionality and use of the temporary pacemaker, based on needs and hospital procedures.

1.5. Environment of use

The DM must be used in hospital environments: cardiology, emergency room and cardiac surgery. The temporary pacemaker must be used in an environment where continuous monitoring of the patient is guaranteed, in order to ensure optimal functioning of the device and guaranteeing that the patient receives the most adequate treatment.

1.6. Essential Performance

The essential performances of the MD can be summarized in two fundamental functions.

- Pacing: the MD must stimulate the heart tissue (with amplitude and frequency that can be set by the user and in the chosen mode - synchronous or asynchronous)
- Sensing: the MD must recognize spontaneous cardiac activity and synchronize with it if the device has been set with ON DEMAND functionality (synchronous).

1.7. Clinical Benefit

Temporary cardiac stimulation is used in the treatment of heart disorders to:

- Heartbeat recovery, keeping the patient alive until spontaneous recovery or until the implantation of a permanent pacemaker
- Maintenance of the correct heart rate in case of pathological alterations of the spontaneous generation of cardiac impulse, or stabilization of the patient from a hemodynamic point of view, i.e. guaranteeing the correct flow of blood into the patient

1.8. Intended purpose

A temporary pacemaker is used in connection to a system of electro-catheters, so to provide a temporary single - chamber stimulation in a hospital environment by trained personnel. The temporary pacemaker can provide either synchronous (on demand) or asynchronous, short-term stimulation for therapeutic, prophylactic or diagnostic purposes, in particular:

- support, management, assessment of a patient with cardiac complications during invasive procedures
- patient support after heart surgery

1.9. Intended patient population

Patients who need short-term, synchronous or asynchronous cardiac electro-stimulation for therapeutic, prophylactic or diagnostic purposes.

2. BATTERY

The device is battery powered. 2 main batteries: AA size 1.5 V alkaline (type LR6). Assessment of the main and stand-by batteries charge is made during the self-test (see the relevant paragraph) and during normal operation.

Pay much attention to the polarity of batteries upon their insertion. The stimulator is protected against damages caused by improper polarity. However, when the replacement is made during the stimulation, the unit would be turned off upon resuming power from the main batteries.

It is recommended to use alkaline batteries (as indicated in § 6) unless the autonomy could be lower than indicated.

Before each use, especially if prolonged, it is firmly recommended replacing the main (alkaline) batteries. The battery replacement will be indispensable if the stimulator has been turned off after the request to change batteries.

It will be necessary to remove the batteries from their compartment if the stimulator is stored for prolonged

periods (longer than 4 weeks).

3. POTENTIAL COMPLICATIONS – RESIDUAL RISKS

Complications likely to occur upon using heart stimulator 4429 Easypace are the following:

- muscle stimulation causing contractions and pain;
- breaking or displacement of the catheter causing intermittence or lack of stimulation
- erroneous setting of stimulation parameters causing capture failure;
- erroneous adjustment of the stimulation rate causing induction and tachycardia;
- erroneous setting of sense parameters causing the stimulation to loose synchronism or be inhibited.

4. CONTRAINDICATIONS

There are no contraindications in relation to the use of 4429 Easypace for temporary cardiac stimulation during a therapy or arrhythmia prevention. However the Patient's state of health or the presence of electromagnetic interference can limit the choice of modes and stimulation parameters.

The stimulator is provided with numerous and redundant safety circuits whose simultaneous failure - unsigned to the User - is an event very unlikely to occur. However, it is recommended that the use of the heart stimulator 4429 Easypace be restricted to authorised medical staff and always accompanied by ECG monitoring.

5. WARNING AND PRECAUTIONS

Indicated here below are the precautions to be taken for the correct use of the heart stimulator 4429 Easypace. Further information can be found in other sections of this manual.

The manufacturer shall hold responsible for the safety, reliability and performance of the unit provided the following indications are strictly observed.

- The device must be used by physicians with extensive knowledge of cardiac stimulation.
- The device must only be used in compliance with the instructions for use. It is necessary to read every part of the manual, if possible in front of the device which is turned on and operating, in order to fully understand its functioning.
- Before each use, especially if prolonged, it is firmly recommended replacing the main (alkaline) batteries. The battery replacement will be indispensable if the stimulator has been turned off after the request to change batteries.
- In the case of prolonged use it will be necessary to make periodic checks of the stimulation parameters (sensing and pacing thresholds).
- It is indispensable to continuously monitor the patient during stimulation. In emergency situations it will also be necessary to have a defibrillator on hand and ready for use.
- The programming of a low value for the sensing threshold increases the probability of the 4429 Easypace operations being influenced by muscular tremors. The device may interpret false signals as spontaneous cardiac activity, with consequent inhibition of the stimulation.
- All the automatic functions offered by the 4429 device have been implemented to help the User make a rapid search for the correct parameters. However, it will always be the physician's responsibility to evaluate the suitability of these values in each specific case.
- In order to use components in contact with the Patient that have to be connected up to the device, follow the instructions supplied with the same.
- In general, an endocavitary or epimiocardiac electrode connected to any type of electro-medical device represents a pathway for direct current with low resistance towards the myocardium. The risk of inducing threatening ventricular arrhythmias resulting from the dispersion of alternate current increases considerably when a device supplied by the mains is connected to a derivation system. The utmost care must therefore be taken in earthing devices supplied from the mains and used in the vicinity of the Patient.
- The heart stimulator 4429 Easypace may also be used within another system consisting of other equipment. In this case the ascertaining of the compliance of the equipment chosen and the specific system with the electrical safety and electromagnetic compatibility provisions in force will be the responsibility of the designer and manufacturer of the same.

In this eventuality it is always recommended taking great care in connecting all the other equipment used in the procedure (electrocardiographs or polygraphs, monitors, etc.) up to the equipotentiality node of the electrical system.

- It is also necessary to avoid touching any metallic parts of the terminals with bare hands and/or letting them come in contact with wet or electrically conductive surfaces.
 - As it is always good practice to equalise the potential between User and Patient, the User is advised to touch the patient in a area that is not too close to the stimulation derivations.
 - Keep all possible sources of static electricity away from the stimulation system.
 - The presence of an electrocatheter in contact with the heart wall may allow for the radiofrequency generated by an electrotome to reach the heart wall itself, producing possible arrhythmias or burns to the myocardium.
 - Devices that generate ionising radiations (cobalt machines, linear accelerators, etc.) could damage the electronic components inside the stimulator, especially following an accumulation of a considerable quantity of radiations. In case of use in the above mentioned conditions it will be necessary to provide an efficacious shield for the stimulator plus continuous monitoring of the Patient.
 - The stimulator is protected from defibrillation discharges at the input. It is however recommended wherever possible disconnecting the electrocatheter from the device before defibrillating.
 - Before each use in heart surgery, check that the locking jaws of the two sockets of the heart surgery cable close evenly and firmly on the terminals. In case of doubt do not use the device and contact FIAB.
 - Before each use check the integrity of the device and all its accessories. Particular attention must be paid when inspecting the cables and connectors. In case of doubt do not use the device and contact FIAB.
 - Before each use check the compatibility of the pin of the electrocatheter with the socket of the device.
- In the event of using electrocatheters not fitted with protected pins, but instead with pins with male contacts with a diameter of less than 2 mm, it is recommended using the heart surgery cable F7818/30 with screw-on sockets in order to allow for stimulation.
- Insert the terminals of the electrode fully into each socket without forcing, following the indications in chapter 10 of IFU - "Connection to the patient and Instruction for use".
 - Before each use check that the transparent protective cover of the keypad is correctly inserted and that there are no breaks or tears. This protection must always be put back into place after setting the parameters.
 - Should the device be dropped, check its correct functioning and in case of doubt contact FIAB before reusing. If there is any apparent mechanical damage avoid using the device and contact FIAB.
 - All operations to open the device, gain access to the internal parts, and repair the same must only ever be carried out by authorized FIAB personnel.
 - The environment where the device is used must be compliant with the safety regulations.
 - The device must never be immersed in water or any other liquids.
 - It will be necessary to remove the batteries from their compartment if the stimulator is stored for prolonged periods (longer than 4 weeks).
 - After storing the 4429 device in environments with temperatures and/or humidity unsuited to its use, before using again wait until a suitable temperature and humidity have been reached as indicated on the "characteristics" card.
 - Product not sterile and not sterilizable: the hygienic conditions of the device cannot be guaranteed in the event of being placed inside a sterile field.
 - In order to avoid damaging the device it must not be sterilized. To clean the device and the accessories follow the relative indications in chapter 12 of IFU.
 - The accessories supplied with the device are not sterile.

IT IS RECOMMENDED HAVING THE 4429 DEVICE CHECKED AT LEAST ONCE A YEAR BY OUR MEASURING LABORATORY.

6. CHECKS

It is advisable a periodic check by the User on the device functions, its accessories and the state of connections. If any failures or bad state of connections are detected, do not use the device but contact FIAB. Moreover it is necessary having the 4429 device checked at least once a year by our measuring laboratory. Our technicians will do the following actions:

- replacing stand-by batteries;
- functioning's check;
- electrical safety test;

- certified final inspection.


7. CLEANING

The User should take care of the cleaning of the unit to meet ordinary sanitary requirements. For the cleaning of the external parts of the unit and accessories, it is recommended using ethyl alcohol or Citrosil®. Belts can be washed with water and soap at a temperature of 40 °C.

For the cleaning of the parts in contact with the Patient, reference is made to the relevant instructions.


8. CONNECTION TO PATIENT

The 4429 Easypace comes with protected outlet sockets, the geometry of which is illustrated in the following figure, and is capable of accepting 2 mm pins, both protected and unprotected. To avoid problems one must be sure that the protected jacks of the electro-catheters in use are compatible with the output bushings, so that they can be inserted easily and guarantee a good electrical contact and mechanical stability.

The protected jack of the electro-catheter or of the extension cable, if any, must possess the following characteristics:

- insulator with internal diameter of 4.5mm and outer diameter of at most 6.2mm (suggested thickness 0.8mm)
- internal contact diameter 2mm

The electrocatheters produced by FIAB that guarantee full compatibility with the 4429 output bushings belong to the "Spike" family.

If the device has to be connected to an epicardial needle it is necessary to use a sterile bipolar adapter cable for heart surgery (e.g. F7818/30) with bushings to tighten.

These bushings can receive sockets with diameter between 0.5mm and 2.0mm. For an ideal coupling, it is necessary that the bushel jaws maximal opening is only slightly larger than the diameter of the socket to be inserted. To regulate the opening, note that it widens turning the bushel external shell counterclockwise and it shrinks in the opposite direction.

The epicardial needle has normally a diameter of about 1mm. Before inserting it in the bushel it is necessary to check that the opening is regulated suitably. In order to ensure the efficacy of the tightening system it is also recommended to insert the needle keeping it oriented as the bushel axis and then tighten the bushel firmly turning it clockwise, to guarantee that the contact is electrically efficient and mechanically stable. Finally, it is recommended to protect the possible residual exposed part of the needle with an insulator.

9. TECHNICAL CHARACTERISTICS

Producer	FIAB SpA Via Costoli, 4- 50039 Vicchio - Firenze (ITALIA)
Model	4429
Class and Type	Portable device, internally powered
Applied part	CF type, protected from defibrillator discharge
Employment conditions	Device for continuous operation
Compliance to safety standards	EN 60601-1:2006 Medical electrical equipment - Part 1: General requirements for basic safety and essential performance EN IEC 60601-2-31:2010 Medical electrical equipment - Part 2-31: Particular requirements for the basic safety and essential performance of external cardiac pacemakers with internal power source EN 60601-1-2:2010 Part 1-2: General requirements for basic safety and essential performance - Collateral Standard: Electromagnetic disturbances - Requirements and tests MDD 93/42/CEE
Intended use	Endocavitary temporary single chamber stimulation
Operation conditions	Temperature: from +10°C to +40 °C Relative humidity: from 30% to 70% RH not condensed Atmospheric pressure: from 700 mbar to 1060 mbar
Storing and transportation conditions (4 weeks max)	Temperature: from 0°C to +50 °C Relative humidity: from 20% to 80% RH not condensed Atmospheric pressure: from 500 mbar to 1060 mbar Packaging suitable to prevent shock induced mechanical damages

Rev.	Revision Status / Description	Date
A	First Issued	12/02/2025

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