NASAL CANNULA FOR OXYGENTHERAPY

INFORMATION FOR SAFETY AND PERFORMANCE

(EN)

INDICATIONS

Therapeutic administration of oxygen through the nose.
Disposable product (for a single patient), to be used continuously for a period not exceeding 30 days.
The oxygen source must be in according to the regulations in force.
The products must be used by qualified healthcare personnel.

CONTRA-INDICATIONS

No contraindications for oxygen therapy are described.

PRECAUTIONS AND WARNINGS

Reuse of the product may result in:

- 1 presence of biological residues that might cause cross-infections,
- 2 alteration of materials,
- 3 loss of initial functional features of product.

Before connecting to the oxygen source, check its operation.

Accidental crushing or kneeling of the connection tubing do not allow oxygen to flow properly and may cause the disjunction of the connector. It is therefore recommended that the users control carefully.
The same attention must be paid even to the models which are supplied with a sure-flow tubing for which the risk of disjunction exists even though in a reduced way.

In case of product containing phthalates (DEHP). Even if risks to human health are not confirmed, the use of the devices on pregnant/nursing or child should be evaluated by your doctor.

The concentrated oxygen can ignite combustible materials. To avoid risk of fire or explosion keep away from possible sources of ignition, and refer to the precautions for use of the equipment used for the delivery of oxygen.

POTENTIAL COMPLICATIONS

The use of nasal cannulas may be associated with epistaxis in subjects prone to bleeding and/or on anticoagulant therapy.

Rev.	Revision Status / Description	Date
A	First Issued	08/05/2023
B	Introduction of form 99501185	03/07/2023

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