

FABT[®]

PRODUCT CATALOGUE

DEFIBRILLATION
+EURODEFIPADS[®]+



EURODEFIPADS®

MULTIFUNCTION PADS



High quality materials, continuous process control and extreme care in assessing the compatibility and the conformity to the international standards make EURODEFIPADS® excellent and safe devices. Easy to be used, they assure optimal clinical performances both in emergency and elective cases.

EURODEFIPADS® is a complete line of adult and pediatric disposable multifunction electrodes for defibrillation, synchronized cardioversion, external cardiac stimulation and ECG monitoring with different connectors compatible with the main defibrillators present in the market.


Radiotransparent models (with radiolucent cable and pads) have been designed for use in EP labs, haemodynamics and during radiographic and angiographic procedures. Lead-out preconnectable models (marked with suffix W) have cable and terminal connector outside the pouch, allowing to connect pads to the defibrillator before opening pouch; they have been intended especially for emergency situations.

Euro Defi Pads are compatible with monophasic and biphasic defibrillators (manual, semi-automatic and automatic) and they comply with **European Directive and International Standards** (IEC/EN 60601-1, IEC/EN 60601-2-4, ANSI/AAMI DF-80, EN ISO 10993-1)

ADULT



PAEDIATRIC



RADIOTRANSSPARENT



LEAD OUT





Please contact the International Sales Dept. in order to get the compatibility list and find out the right reference for each defibrillator brand and model. Brands, models and commercial names reported in this document belong to the companies owning the corresponding registered brands, models and names. FIAB is not affiliated to any of the above mentioned companies.

CONNECTORS BRANDS

- F7950 F7950P F7950PL
- F7950W F7950PW F7950WPL
- F7750 F7750P
- F7750W
























Philips Medical®
Laerdal Medical®
Agilent Philips Medical













ADULT

UNIVERSAL

PAEDIATRIC

	Radiopaque	Radiotransparent		Radiopaque	AED	Radiotransparent
						
						
F7950 F7950W	F7750 F7750W	F7950U F7950UW F7750U	F7950P F7950PW	F7950PL F7950WPL	F7750P	

						
						
Zoll Medical Corp® M&B	F7951 F7951W F7951/300	F7751 F7751W		F7951P F7951PW		F7751P F7751PW

							
							
Physiocontrol® Bexen Cardio Cardioline® Mindray Comen	F7952 F7952W F7952/200	F7752 F7752W		F7952P F7952PW	F7952PL F7952WPL	F7752P F7752PW	

EURODEFIPADS®

Please contact the International Sales Dept. in order to get the compatibility list and find out the right reference for each defibrillator brand and model.

Brands, models and commercial names reported in this document belong to the companies owning the corresponding registered brands, models and names. FIAB is not affiliated to any of the above mentioned companies.

ADULT

UNIVERSAL

PAEDIATRIC

Radiopaque

Radiotransparent

Radiopaque

AED

Radiotrasparent



CONNECTORS	BRANDS	ADULT Radiopaque	ADULT Radiotransparent	UNIVERSAL	PAEDIATRIC Radiopaque	PAEDIATRIC AED	PAEDIATRIC Radiotrasparent
	Agilent Laerdal Medical® Hewlett-Packard	F7953 F7953W	F7753 F7753W		F7953P F7953PW		F7753P
	Schiller Esaote	F7954 F7954W	F7754		F7954P		F7754P
	Nihon Kohden Welch Allyn® Innomed MRL	F7955 F7955W	F7755 F7755W		F7955P F7955PW		F7755P
	F7956 F7956W F7756 F7956P F7956PW F7756P Schiller	F7956 F7956W	F7756		F7956P F7956PW		F7756P
	Cardiac Science	F7957W				F7957PL	
	CU Medical Systems, inc	F7958 F7958W			F7958P F7958PW	F7958PL	
	CU Medical Systems, inc Cmos Drake	F7959 F7959W	F7759		F7959P	F7959PL F7959WPL	
	Metrax Primedic	F7960W	F7760W				
	Nihon Kohden		F7761	F7961U F7961UW F7961UWXL F7761U	F7961P F7961PW		
	Marquette GE® Esaote	F7962	F7762		F7962P		
	Nihon Kohden	F7963					
	Weinmann Cardiaid	F7964W			F7964PW		

Please contact the International Sales Dept. in order to get the compatibility list and find out the right reference for each defibrillator brand and model.

Brands, models and commercial names reported in this document belong to the companies owning the corresponding registered brands, models and names. FIAB is not affiliated to any of the above mentioned companies.

		ADULT		UNIVERSAL	PAEDIATRIC		
		Radiopaque	Radiotransparent		Radiopaque	AED	Radiotransparent
CONNECTORS	BRANDS						
	Corpuls	F7965 F7965/PLUS			F7965P		
	Corpuls	F7965W/NG F7965/PLUSNG			F7965PW/NG		
	Defibtech	F7966W				F7966WPL	
	Schiller	F7967 F7967W			F7967P F7967PW		
	Mediana	F7969W			F7969PW		
	GE®	F7981			F7981P		
	Mediana Tecnogaz	F7986 F7986W		F7986UW F7986UW/D700 F7786UW/D700	F7986P F7986PW		
	Zoll Medical Corp®	F7987W			F7987PW		
	Nihon Kohden Not for sale in France, Germany and UK			F7989UW			
	Radian			F7992U			
	Laerdal Medical®	F7970			F7970P		
	Physiocontrol®	F7971					
	Marquette GE® Esaote	F7972			F7972P		




EURODEFIPADS®

STIMULATION PADS

CE 2797

Besides EURODEFIPADS® models, some additional pads have been developed only for external cardiac stimulation.



CONNECTORS	BRANDS		
	Zoll Medical Corp® (M series, pacer version only)	adult paediatric	F7951/STIM F7951P/STIM
	Physiocontrol®	adult	F7973
	Nihon Kohden Innomed	adult	F7955/STIM


TRAINER PADS

CE

Defibrillation training pads (adult and paediatric) are equipped with snap connection and must be used with the appropriate connection cable. The pad's special shape makes the application extremely easy both on manikins and on other surfaces.

PADS

F7999	adult
F7999P	paediatric



CONNECTING CABLES

CONNECTORS	BRANDS	
	Philips Medical® Laerdal Medical® Agilent Philips Medical	F5350
	Physiocontrol® Osatu Bexen Cardio Cardioline® - Mindray	F5352
	Welch Allyn	F5355
	Schiller	F5356
	Weinmann Cardiaid	F5364
	Defibtech	F5366
	Mediana	F5369
	Mediana Tecnogaz	F5386
	Nihon Kohden	F5389



FIAB SpA
Via Costoli, 4 - I-50039 Vicchio (Florence) - Italy
Tel. +39 0558497999 - Fax +39 055844381
export@fiab.it - www.fiab.it

