

**FLAB**<sup>®</sup>

PRODUCT CATALOGUE

# DEFIBRILLATION **+EURODEFIPADS<sup>®</sup>**





**EURODEFIPADS®**

## MULTIFUNCTION PADS



*Manufactured with high-quality materials, subjected to strict and continuous process control and thoroughly tested for compatibility and compliance with international standards, EURODEFIPADS® ensure both safety and excellence.*

*Easy to use, they deliver optimal clinical performance in both emergency situations and elective procedures.*

EURODEFIPADS® is a complete line of adult and pediatric disposable multifunction electrodes for defibrillation, synchronized cardioversion, external cardiac stimulation and ECG monitoring with different connectors compatible with the main defibrillators present in the market.

Radiotransparent models (with radiolucent cable and pads) have been designed for use in EP labs, haemodynamics and during radiographic and angiographic procedures. Lead-out preconnected models (marked with the suffix W) have the cable and connector placed outside the pouch, so the pads can be connected to the defibrillator before opening the packaging. This design is especially useful in emergency situations.

EURODEFIPADS® are classified as Class IIb medical devices in accordance with the Medical Device Regulation (EU) 2017/745 (MDR) and are CE marked following assessment by a Notified Body, ensuring full compliance with the highest standards of safety, performance and quality.

### ADULT



### PAEDIATRIC



### RADIOTRANSSPARENT



### LEAD OUT



# EURODEFIPADS®



Please contact the International Sales Dept. in order to get the compatibility list and find out the right reference for each defibrillator brand and model.

Brands, models and commercial names reported in this document belong to the companies owning the corresponding registered brands, models and names. FIAB is not affiliated to any of the above mentioned companies.

## ADULT

## UNIVERSAL

## PAEDIATRIC

ADULT		UNIVERSAL	PAEDIATRIC		
Radiopaque	Radiotransparent		Radiopaque	AED	Radiotransparent

### CONNECTORS BRANDS
















CONNECTORS	BRANDS	ADULT Radiopaque	ADULT Radiotransparent	UNIVERSAL	PAEDIATRIC Radiopaque	PAEDIATRIC AED	PAEDIATRIC Radiotransparent
F7950P F7950PW F7750P 	Philips Medical® others	F7950 F7950W	F7750 F7750W	F7950U F7950UW F7750U	F7950P F7950PW	F7950PL F7950WPL	F7750P
	Zoll Medical Corp® M&B	F7951 F7951W F7951/300	F7751 F7751W		F7951P F7951PW		F7751P F7751PW
F7952PL F7952WPL 	Physiocontrol® Bexen Cardio Cardioline® Mindray Comen Jousing Acoresmed others	F7952 F7952W F7952/200	F7752 F7752W		F7952P F7952PW	F7952PL F7952WPL	F7752P F7752PW
	Agilent Laerdal Medical® Hewlett-Packard	F7953 F7953W	F7753 F7753W		F7953P F7953PW		F7753P
	Schiller Esaote	F7954 F7954W	F7754		F7954P		F7754P
	Nihon Kohden Welch Allyn® Artema Innomed MRL Artema S&W Dräger	F7955 F7955W	F7755 F7755W		F7955P F7955PW		F7755P

# EURODEFIPADS®

Please contact the International Sales Dept. in order to get the compatibility list and find out the right reference for each defibrillator brand and model.

Brands, models and commercial names reported in this document belong to the companies owning the corresponding registered brands, models and names. FIAB is not affiliated to any of the above mentioned companies.

**CONNECTORS**      **BRANDS**

		ADULT		UNIVERSAL	PAEDIATRIC		
		Radiopaque	Radiotransparent		Radiopaque	AED	Radiotrasparent
 <p>F7956 F7956W F7756</p> <p>F7956P F7956PW F7756P</p> <p>Schiller</p>		 <p>F7956 F7956W</p>	 <p>F7756</p>		 <p>F7956P F7956PW</p>		 <p>F7756P</p>
 <p>Cardiac Science GE®</p>		F7957W				F7957PL	
 <p>CU Medical Systems, inc</p>		F7958 F7958W			F7958P F7958PW	F7958PL	
 <p>CU Medical Systems, inc Cmos Drake Amoul® Ambulanc (Shenzhen) TECH. CO., LTD</p>		F7959 F7959W	F7759		F7959P	F7959PL F7959WPL	
 <p>Metrax Primedic</p>		F7960W	F7760W				
 <p>Nihon Kohden</p>			F7761	F7961U F7961UW F7961UWXL F7761U	F7961P F7961PW		
 <p>Marquette GE® Esaote</p>		F7962	F7762		F7962P		
 <p>Nihon Kohden</p>		F7963					
 <p>Weinmann</p>		F7964W			F7964PW		
 <p>Corpuls</p>		F7965			F7965P		
 <p>F7965W/NG F7965/PLUSNG</p> <p>F7965PW/NG</p> <p>Corpuls</p>		F7965W/NG F7965/PLUSNG			F7965PW/NG		

# EURODEFIPADS®

Please contact the International Sales Dept. in order to get the compatibility list and find out the right reference for each defibrillator brand and model.

Brands, models and commercial names reported in this document belong to the companies owning the corresponding registered brands, models and names. FIAB is not affiliated to any of the above mentioned companies.

ADULT

UNIVERSAL

PAEDIATRIC

Radiopaque

Radiotransparent

Radiopaque

AED

Radiotrasparent




CONNECTORS

BRANDS

CONNECTORS	BRANDS	Radiopaque	Radiotransparent	UNIVERSAL	Radiopaque	AED	Radiotrasparent
	Defibtech	F7966W				F7966WPL	
	Schiller Philips	F7967 F7967W			F7967P F7967PW		
	Mediana	F7969W			F7969PW		
	GE®	F7981			F7981P		
	Mediana Tecnogaz	F7986 F7986W		F7986UW F7986UW/D700 F7786UW/D700	F7986P F7986PW		
	Zoll Medical Corp®	F7987W			F7987PW		
	Nihon Kohden Not for sale in France, Germany and UK			F7989UW			
	Radian			F7992U			
	Laerdal Medical®	F7970			F7970P		
	Physiocontrol®	F7971					
	Marquette GE® Esaote	F7972			F7972P		

Besides EURODEFIPADS® models, some additional models are intended for external cardiac stimulation only.



CONNECTORS	BRANDS		
	Zoll Medical Corp® (M series, pacer version only)	adult paediatric	F7951/STIM F7951P/STIM
	Physiocontrol®	adult	F7973
	Nihon Kohden Innomed	adult	F7955/STIM

## TRAINER PADS

Defibrillation training pads (adult and paediatric) are equipped with snap connection and must be used with the appropriate connection cable. The pad's special shape makes the application extremely easy both on manikins and on other surfaces.

### PADS

REF.	PATIENT
F7999	adult
F7999P	paediatric



### CONNECTING CABLES

CONNECTORS	BRANDS	
	Philips Medical® Laerdal Medical® Agilent Philips Medical	F5350
	Physiocontrol® Osatu Bexen Cardio Cardioline® - Mindray	F5352
	Welch Allyn	F5355
	Schiller	F5356
	Weinmann Cardiaid	F5364
	Defibtech	F5366
	Mediana	F5369
	Mediana Tecnogaz	F5386
	Nihon Kohden	F5389



FIAB SpA  
Via Costoli, 4 - I-50039 Vicchio (Florence) - Italy  
Tel. +39 0558497999  
export@fiab.it - www.fiab.it



MD 77846