



PRODUCT CATALOGUE

ECG

- DISPOSABLE ELECTRODES
- REUSABLE ELECTRODES
- PATIENT CABLES
- ADAPTERS
- HOLTER CABLES
- GEL

















DISPOSABLE ELECTRODES

EUROEGCELECTRODES™

FIAB EUROEGCELECTRODES provide the highest quality and performance during ECG test, stress test, long and short monitoring, holter test. FIAB electrodes are available in different shapes and sizes. Solid gel has a special formula that improves adhesion and easy application. Liquid gel electrodes allow a quick and stable ECG tracing. Adhesive backing is made in foam or not woven. Connection are standard snap or 4 mm banana adapter. Carbon snap electrodes are excellent for use with RX and MRI. A complete line of prewired electrodes are specially designed for sensitive skin of neonates and paediatric patients. Conductive adhesive provides an excellent ECG reading and the soft medical foam or the not woven with a sensitive adhesive ensure a secure comfortable contact.

SNAP/SOCKET ELECTRODES

ADULT

ref.	picture	suggested use	ref.	picture	suggested use
F9060 48x50mm FOAM solid gel snap		ECG recording long term monitoring holter stress test resting ECG	F9049/100 34x51mm FOAM solid gel snap		ECG recording long term monitoring holter stress test resting ECG
F9079/100 36x40mm FOAM solid gel snap		ECG recording holter stress test resting ECG	F9069RM 42x45mm FOAM solid gel gel snap radiotrasparent		ECG recording resting ECG MRI
F9089/100 36x45mm FOAM solid gel snap		ECG recording holter stress test resting ECG	F9046 28x44mm FOAM solid gel gel snap		ECG recording resting ECG
F2080 Ø 55mm NON WOVEN solid gel snap		ECG recording holter long term monitoring stress test resting ECG	F9046RM 28x44mm FOAM solid gel gel snap radiotrasparent		ECG recording resting ECG MRI
F2017 Ø 50 mm FOAM liquid gel snap		ECG recording holter long term monitoring stress test resting ECG	F9047 32x44mm FOAM liquid gel snap		ECG recording holter stress test resting ECG
F9047/4F 34x57mm FOAM solid gel 4mm socket		ECG recording holter stress test resting ECG	F9047/4FL 28x44mm FOAM liquid gel 4mm socket		ECG recording holter stress tes resting ECG
F9048 34x57mm FOAM solid gel snap		ECG recording holter stress test resting ECG	F9067 48x50mm FOAM liquid gel snap		ECG recording holter stress test resting ECG

PAEDIATRIC

ref.	picture	suggested use	ref.	picture	suggested use
F9060P Ø 40mm FOAM solid gel snap		ECG recording long term monitoring holter stress test resting ECG	F9089P/100 32X36mm FOAM solid gel snap		ECG recording holter stress test resting ECG
F2080P Ø 40mm NON WOVEN solid gel snap		ECG recording long term monitoring holter stress test resting ECG	PG10C Ø 26mm FOAM solid gel snap		ECG recording resting ECG
F2017P Ø 42mm FOAM liquid gel snap		ECG recording long term monitoring holter stress test resting ECG	F9069PRM 32X36mm oval FOAM solid gel snap radiotrasparent		ECG recording resting ECG MRI
F9079P/100 32X36mm FOAM solid gel snap		ECG recording holter stress test resting ECG	F9069NRM Ø 26mm FOAM solid gel snap radiotrasparent		ECG recording resting ECG MRI

TAB ELECTRODES



ADULT

PAEDIATRIC

ref.	picture	suggested use
F3001ECG 23X34 mm PLASTIC solid radiotrasparent		resting ECG

ref.	picture	suggested use
F3001NECG 13X34 mm PLASTIC solid radiotrasparent		resting ECG

PREWIRED ELECTRODES

ref	dimensions	backing	lead length	connector	gel	radioltransparent	picture		
FOR ADULTS									
F9053	Ø 32 mm	FOAM	35 cm	4 mm socket	solid		 <p>F9053</p>		
F9053/10			10 cm					NO	
F9053RM			35 cm					YES	
F9053RM/80			80 cm					YES	
PAEDIATRIC									
F9058N	23 X 30 mm	FOAM	50 cm	1,5 mm DIN socket	solid		 <p>F9053N</p>		
F9058N/F2	23 X 30 mm		50 cm	2 mm socket				NO	
F9058N/60	23 X 30 mm		60 cm	1,5 mm DIN socket				NO	
F9058NRM	23 X 30 mm		50 cm					NO	
F9058/90NRM	23 X 30 mm		90 cm					YES	
F9058P	Ø 32 mm		50 cm					NO	
F9058P/60	Ø 32 mm		60 cm					NO	
F9058PR	20 X 15 mm		50 cm					NO	
F9058PRM	Ø 32 mm		50 cm					YES	
F9058PRRM	20 X 15 mm		50 cm					YES	
F9053N	23 X 30 mm		35 cm					4 mm socket	NO
F9053N/10	23 X 30 mm		10 cm						NO
F9053PR	20 X 15 mm		35 cm	NO					
F9053PR/10	20 X 15 mm		10 cm	NO					
F9053NRM	23 X 30 mm		35 cm	YES					
F9053NRM/80	23 X 30 mm		80 cm	YES					
F9053PRRM	20 X 15 mm	10 cm	YES						
F9054N	23 X 30 mm	35 cm	NON WOVEN	1,5 mm DIN socket	NO				
F9054N3/10	23 X 30 mm	10 cm			NO				
F9054N3/55	23 X 30 mm	55 cm			NO				
F9054NQ/RM/80	22 X 22 mm	80 cm			YES				
F9057P	Ø 32 mm	50 cm			NO				
F9057PRM	Ø 32 mm	50 cm			YES				
F9057N	23 X 30 mm	50 cm			NO				
F9057NRM	23 X 30 mm	50 cm			YES				
F9057PR	20 X 15 mm	50 cm			NO				
F9057PRRM	20 X 15 mm	50 cm			YES				
F9057PQ/60	32 X 32 mm	60 cm			NO				
F9057NQ/60	22 X 22 mm	60 cm			NO				
F9057NQ/RM/50	22 X 22 mm	50 cm	YES						
 <p>F9058P/60</p>  <p>F9058NRM</p>  <p>F9058</p>  <p>F9054</p>  <p>F9057NQ/60</p>									

REUSABLE ELECTRODES

The range of reusable electrodes includes peripheral clamps, peripheral plates with rubber belts, suction chests and spare parts.

PERIPHERAL ELECTRODES

ref	description	picture
CLAMPS FOR ADULTS		
F9024SSC	clamp electrode Silver/Silver Chloride	
F9024NI	clamp electrode Nickel plated	
F9024OSSC	clamp electrode Silver/Silver Chloride, for Obese	
F9024ONI	clamp electrode Nickel plated, for obese	
PAEDIATRIC CLAMPS		
F9023SSC	clamp electrode Silver/Silver Chloride	
F9023NI	clamp electrode Nickel plated	
PLATES		
F9010SSC	Silver/Silver Chloride for adults	
F9010PSSC	Silver/Silver Chloride for paediatric	
F9011PG	rubber belt for peripheral plate electrode h=24mm, L=48cm, for adults	
F9011P	rubber belt for peripheral plate electrode h=17mm, L=26cm, for paediatric	

SUCTION CHEST ELECTRODES


ref	description	picture
WITH RUBBER RING		
F9002SSC	Ø 24 mm, Silver/Silver Chloride	
F9002NI	Ø 24 mm, Nickel plated	
F9003SSC	Ø 30 mm, Silver/Silver Chloride	
F9003NI	Ø 30 mm, Nickel plated	
WITH METALLIC RING		
F9008SSC	Ø 15 mm, Silver/Silver Chloride	
F9008NI	Ø 15 mm, Nickel plated	
F9015SSC	Ø 20 mm, Silver/Silver Chloride	
F9009SSC	Ø 24 mm, Silver/Silver Chloride	
F9009NI	Ø 24 mm, Nickel plated	
F9016SSC	Ø 30 mm, Silver/Silver Chloride	
F9016NI	Ø 30 mm, Nickel plated	
SPARE PARTS		
F9005SSC	spare metallic chest, Ø 15mm, Silver/Silver Chloride	
F9005/20SSC	spare metallic chest, Ø 20mm, Silver/Silver Chloride	
F9005/24SSC	spare metallic chest, Ø 24mm, Silver/Silver Chloride	
F9006SSC	spare metallic chest, Ø 30mm, Silver/Silver Chloride	
F9007	spare rubber ball chest for series: F9003, F9009, F9016	
F9007P	spare rubber ball chest for series: F9002, F9008, F9015	

PATIENT CABLES

FIAB manufactures a wide range of ECG patient cables with direct connection to the main part of devices present on the market. Patient cables are provided with 4 mm plugs, snaps or pinches.



CABLES WITH MOULDED MANIFOLD






CABLES WITH INSPECTABLE MANIFOLD



ECG DEVICE MANUFACTURER

	ASPEL BIONET BTL CARDIOLINE CARDIOTEK CARDIORAPID DR.LEE		BURDICK / QUINTON		ESAOTE
	BOSCH CAMBRIDGE KENZ NIHON KOHDEN PICKER		CARDIOLINE		FUKUDA DENSHI
	BOSCH CUSTO MED GMBH HELLIGE SIEMENS		CARDIOLINE KENZ		FUKUDA DENSHI
	BURDICK / QUINTON SIEMENS		CARDIOLINE CARDIETTE		HELLIGE MARQUETTE
	BURDICK / QUINTON SIEMENS		ESAOTE		MORTARA

SPARE PARTS

ref.	description	picture
KITS OF 10 UNIPOLAR LEADS		
F5377/KIT	4mm plug connectors, without protection	
F5388/KIT	4mm plug connectors, 4,7 Kohm protection	
F5388R/KIT	4mm plug connectors, 10 Kohm protection	
F5376/KIT	Pinch connectors, without protection	
F5386/KIT	Pinch connectors, 4,7 Kohm protection	
F5386R/KIT	Pinch connectors, 10 Kohm protection	
F5399/KIT	Snap connectors, without protection	

ADAPTERS

A large range of adapters from banana plug to snap or alligator allows to modify the configuration of each patient cables or single leads.

ref.	description	picture
F9017W	From 4mm socket to snap+alligator, white	
F9017/4PW	from 4mm socket to pich, white	
F9019	from 4mm socket to alligator (*) (*) suffix colour key: R=red, N=black	
PG922/4T	from 4mm socket to snap (**) (**) su ffix colour key: R=red, N=black, V=green, G=yellow, B=white	
PG922/4KIT	from 4mm socket to snap , kit of 10 pieces (***) (***) 10pcs: 6 white, 1 black, 1green, 1 yellow	

HOLTER CABLES

Holter cables have a snap connection at one side and different connections at the other one, in order to match with every needing. Several colours are available for every single cable.

ref.	description	picture
SINGLE CABLES		
F9022S	L=60cm, 2mm safety plug (*)	
F9052S	L=80cm, 2mm safety plug (*)	
F9032S	L=150cm, 2mm safety plug (*)	
F9022D	L=60cm, safety 1,5mm DIN socket (*)	
F9052D	L=90 cm, safety 1,5 mm DIN socket (*)	
F9032D	L=120 cm, safety 1,5 mm DIN socket (*)	
F9032DN/200	L=200 cm, safety 1,5 mm DIN socket, black	
CABLES IN KIT		
F9022D/KITX	Kit of dnap connection cables, L=60cm, safety 1,5mm DIN socket (**)	
F9052D/KITX	Kit of dnap connection cables, L=90cm, safety 1,5mm DIN socket (**)	
F9032D/KITX	Kit of dnap connection cables, L=120cm, safety 1,5mm DIN socket (**)	

(*) suffix colour key: R=red, N=black, A= Orange, V=green, G=yellow, W=white, M=brown, B=blue, I=purple

(**) F9022D/KITX and F9052D/KITX models are available with 3, 5, 7 or 8 cables, please replace "X" with the related number

(***) F9032D/KITX is available with 3, 5 or 7 cables, please replace "X" with the related number

GEL

GEL range includes models specifically intended for ECG and ultrasound: dispensers are available separately as spare parts.

ECG

ref.	description	picture
G004	spray gel, 230ml bottle	
G005	liquid gel, 260ml bottle	

ULTRASOUND

ref.	description	picture
G006	liquid gel, 260ml bottle	
G007	liquid gel, 1000ml bottle	
G0064	liquid gel, 5000ml flexible bag, 260ml refilling bottle	
G0066	liquid gel, 5000ml jerry can, 260ml refilling bottle	

ULTRASOUND AND IPL

ref.	description	picture
G002	transparent liquid gel, 20ml pouch	
G008	transparent liquid gel, 260ml bottle	
G009	transparent liquid gel, 1000ml bottle	
G010	transparent liquid gel, high density, 1000ml jar	
G080	transparent liquid gel, high density, 500 ml bottle	
G0084	transparent liquid gel, 5000ml flexible bag, 260ml refilling bottle	

DISPENSERS

ref.	description	picture
G0067	reusable dispenser for 5000ml jerry can	
G0068	reusable dispenser for 5000ml flexible bag	
G0069	reusable dispenser for 1000ml bottle	



FIAB SpA
Via Costoli, 4 - I-50039 Vicchio (Florence) - Italy
Tel. +39 0558497999 - Fax +39 055844381
export@fiab.it - www.fiab.it



MD 77846