



PRODUCT CATALOGUE

ECG

- DISPOSABLE ELECTRODES
- REUSABLE ELECTRODES
- PATIENT CABLES
- ADAPTERS
- HOLTER CABLES
- GEL



DISPOSABLE ELECTRODES













EUROECGELECTRODES™

FIAB EUROECGELECTRODES are designed to deliver maximum quality and performance across a wide range of applications, including ECG tests, stress tests, short- and long-term monitoring and Holter exams.

Available in various shapes and sizes, FIAB electrodes meet the diverse clinical needs of healthcare professionals. The solid gel formulation ensures excellent adhesion and easy application. Liquid gel versions allow for quick, stable and accurate ECG signal acquisition. The adhesive backing is offered in foam or non-woven material, providing flexibility and patient comfort. Electrodes come with standard snap connectors or 4 mm banana adapters for universal compatibility. Carbon snap electrodes are ideal for use during X-ray and MRI procedures, minimizing image artifacts.



FIAB also offers a complete line of prewired electrodes, specially developed for the sensitive skin of neonates and pediatric patients.

Conductive adhesive provides an excellent ECG signal quality and gentle materials (soft medical foam or sensitive non-woven adhesive) ensure secure and comfortable contact.








Electrode type	REF and dimensions	Picture	Suggested Use	REF and dimensions	Picture	Suggested Use
FOAM WITH SOLID GEL	F9060 48x50mm		ECG recording, long term monitoring, holter, stress test, resting ECG	F9079/100 36x40mm		ECG recording, holter, stress test, resting ECG
	F9089/100 36x40mm		ECG recording, holter, stress test, resting ECG	F9046 28x44mm		ECG recording, resting ECG
	F9048 34x57mm		ECG recording, holter, stress test, resting ECG	F9049/100 34x51mm		ECG recording, long term monitoring, holter, stress test, resting ECG
	F9069RM 42x45mm Radio-transparent		ECG recording, resting ECG, MRI	F9046RM 28x44mm Radio-transparent		ECG recording, resting ECG, MRI
FOAM WITH LIQUID GEL	F9067 48x50mm		ECG recording, holter, stress test, resting ECG	F9047 32x44mm		ECG recording, holter, stress test, resting ECG
	F2017 Ø 50mm		ECG recording, long term monitoring, holter, stress test, resting ECG			
NON WOVEN WITH SOLID GEL	F2080 Ø 50mm		ECG recording, long term monitoring, holter, stress test, resting ECG			

EUROECGELECTRODES™


4MM SOCKET ELECTRODES FOR ADULTS


Electrode type	REF and dimensions	Picture	Suggested Use
FOAM WITH SOLID GEL	F9047/4F 28x44mm		ECG recording, holter, stress test, resting ECG
FOAM WITH LIQUID GEL	F9047/4FL 28x44mm		ECG recording, holter, stress test, resting ECG

PAEDIATRIC SNAP ELECTRODES

Electrode type	REF and dimensions	Picture	Suggested Use	REF and dimensions	Picture	Suggested Use
FOAM WITH SOLID GEL	F9060P Ø 40mm		ECG recording, long term monitoring, holter, stress test, resting ECG	F9079P/100 32x36mm		ECG recording, holter, stress test, resting ECG
	F9089P/100 32x36mm		ECG recording, holter, stress test, resting ECG	PG10C Ø 26mm		ECG recording, resting ECG
	F9069PRM 32x36mm Radio-transparent		ECG recording, resting ECG, MRI	F9069NRM Ø 26mm Radio-transparent		ECG recording, resting ECG, MRI
NON WOVEN WITH SOLID GEL	F2080P Ø 40mm		ECG recording, long term monitoring, holter, stress test, resting ECG			






TAB ELECTRODES FOR RESTING ECG

ADULT	
ref.	picture
F3001ECG 23x34mm PLASTIC WITH SOLID GEL RADIOTRANSPARENT	

PAEDIATRIC	
ref.	picture
F3001NECG 13x34mm PLASTIC WITH SOLID GEL RADIOTRANSPARENT	

EUROECGELECTRODES™




PREWIRED ELECTRODES WITH SOLID GEL

FOR ADULTS									
ref.	dimensions	backing	lead length	connector	radiopaque	picture			
F9053	Ø 32 mm	FOAM	35 cm	4 mm socket	NO				
F9053/10			10 cm		NO				
F9053RM			35 cm		YES				
F9053RM/80			80 cm		YES				
PAEDIATRIC									
ref.	dimensions	backing	lead length	connector	radiopaque	picture			
F9058N	23 X 30 mm	FOAM	50 cm	1,5 mm DIN socket	NO				
F9058NRM	23 X 30 mm		50 cm		NO				
F9058N/60	23 X 30 mm		60 cm		NO				
F9058/90NRM	23 X 30 mm		90 cm		YES				
F9058P	Ø 32 mm		50 cm		NO				
F9058P/60	Ø 32 mm		60 cm		NO				
F9058PR	20 X 15 mm		50 cm		NO				
F9058PRM	Ø 32 mm		50 cm		YES				
F9058PRRM	20 X 15 mm		50 cm		YES				
F9053N	23 X 30 mm		FOAM		35 cm		4 mm socket	NO	
F9053N/10	23 X 30 mm				10 cm			NO	
F9053PR	20 X 15 mm				35 cm			NO	
F9053PR/10	20 X 15 mm				10 cm			NO	
F9053NRM	23 X 30 mm				35 cm			YES	
F9053NRM/80	23 X 30 mm				80 cm			YES	
F9053PRRM	20 X 15 mm				10 cm			YES	
F9054N	23 X 30 mm	NON WOVEN		35 cm	1,5 mm DIN socket	NO			
F9054N3/10	23 X 30 mm			10 cm		NO			
F9054N3/55	23 X 30 mm			55 cm		NO			
F9054NQ/RM/80	22 X 22 mm		80 cm	YES					
F9057P	Ø 32 mm		50 cm	NO					
F9057PRM	Ø 32 mm		50 cm	YES					
F9057N	23 X 30 mm		50 cm	NO					
F9057NRM	23 X 30 mm		50 cm	YES					
F9057PR	20 X 15 mm		50 cm	NO					
F9057PRRM	20 X 15 mm		50 cm	YES					
F9057PQ/60	32 X 32 mm	60 cm	NO						
F9057NQ/60	22 X 22 mm	60 cm	NO						
F9057NQ/RM/50	22 X 22 mm	50 cm	YES						

REUSABLE ELECTRODES

The range of reusable electrodes includes peripheral clamps, peripheral plates with rubber belts, suction chests and spare parts.

PERIPHERAL ELECTRODES

CLAMPS FOR ADULTS		
ref.	description	picture
F9024SSC	clamp electrode Silver/Silver Chloride	
F9024NI	clamp electrode Nickel plated	
F9024OSSC	clamp electrode Silver/Silver Chloride, for Obese	
F9024ONI	clamp electrode Nickel plated, for obese	
PAEDIATRIC CLAMPS		
ref.	description	picture
F9023SSC	clamp electrode Silver/Silver Chloride	
F9023NI	clamp electrode Nickel plated	
PLATES		
ref.	description	picture
F9010SSC	Silver/Silver Chloride for adults	
F9010PSSC	Silver/Silver Chloride for paediatric	
F9011PG	rubber belt for peripheral plate electrode h=24mm, L=48cm, for adults	
F9011P	rubber belt for peripheral plate electrode h=17mm, L=26cm, for paediatric	

SUCTION CHEST ELECTRODES

WITH RUBBER RING		
ref.	description	picture
F9002SSC	Ø 24 mm, Silver/Silver Chloride	
F9002NI	Ø 24 mm, Nickel plated	
F9003SSC	Ø 30 mm, Silver/Silver Chloride	
F9003NI	Ø 30 mm, Nickel plated	
WITH METALLIC RING		
ref.	description	picture
F9008SSC	Ø 15 mm, Silver/Silver Chloride	
F9008NI	Ø 15 mm, Nickel plated	
F9015SSC	Ø 20 mm, Silver/Silver Chloride	
F9009SSC	Ø 24 mm, Silver/Silver Chloride	
F9009NI	Ø 24 mm, Nickel plated	
F9016SSC	Ø 30 mm, Silver/Silver Chloride	
F9016NI	Ø 30 mm, Nickel plated	
SPARE PARTS		
ref.	description	picture
F9005SSC	spare metallic chest, Ø 15mm, Silver/Silver Chloride	
F9005/20SSC	spare metallic chest, Ø 20mm, Silver/Silver Chloride	
F9005/24SSC	spare metallic chest, Ø 24mm, Silver/Silver Chloride	
F9006SSC	spare metallic chest, Ø 30mm, Silver/Silver Chloride	
F9007	spare rubber ball chest for series: F9003, F9009, F9016	
F9007P	spare rubber ball chest for series: F9002, F9008, F9015	

PATIENT CABLES

FIAB manufactures a wide range of ECG patient cables with direct connection to the main part of devices present on the market. Patient cables are provided with 4 mm plugs, snaps or pinches.



CABLES WITH MOULDED MANIFOLD

F8000 series
small moulded manifold
4mm plug connection



F8000 series
small moulded manifold
snap connection



CABLES WITH INSPECTABLE MANIFOLD

F6000 series
small manifold
4mm plug connection



F6000P series
small manifold
pinch connection



F5000c series
large manifold
snap connection






ECG DEVICE MANUFACTURER

	ASPEL BIONET BTL CARDIOLINE CARDIOTEK CARDIORAPID DR.LEE		BURDICK / QUINTON		ESAOTE
	BOSCH CAMBRIDGE KENZ NIHON KOHDEN PICKER		CARDIOLINE		FUKUDA DENSHI
	BOSCH CUSTO MED GMBH HELLIGE SIEMENS		CARDIOLINE KENZ		FUKUDA DENSHI
	BURDICK / QUINTON SIEMENS		CARDIOLINE CARDIETTE		HELLIGE MARQUETTE
	BURDICK / QUINTON SIEMENS		ESAOTE		MORTARA





SPARE PARTS

KITS OF 10 UNIPOLAR LEADS

ref.	description	picture
F5377/KIT	4mm plug connectors, without protection	
F5388/KIT	4mm plug connectors, 4,7 Kohm protection	
F5388R/KIT	4mm plug connectors, 10 Kohm protection	
F5376/KIT	Pinch connectors, without protection	
F5386/KIT	Pinch connectors, 4,7 Kohm protection	
F5386R/KIT	Pinch connectors, 10 Kohm protection	
F5399/KIT	Snap connectors, without protection	




ADAPTERS

A large range of adapters from banana plug to snap or alligator allows to modify the configuration of each patient cables or single leads.

ref.	description	picture
F9017W	From 4mm socket to snap+alligator, white	
F9017/4PW	from 4mm socket to pich, white	
F9019	from 4mm socket to alligator (*) (*) suffix colour key: R=red, N=black	
PG922/4T	from 4mm socket to snap (**) (**) suffix colour key: R=red, N=black, V=green, G=yellow, B=white	
PG922/4KIT	from 4mm socket to snap , kit of 10 pieces (***) (***) 10pcs: 6 white, 1 black, 1green, 1 yellow	

HOLTER CABLES

Holter cables have a snap connection at one side and different connections at the other one, in order to match with every needing. Several colours are available for every single cable.

SINGLE CABLES		
ref.	description	picture
F9022S	L=60cm, 2mm safety plug (*)	
F9052S	L=80cm, 2mm safety plug (*)	
F9032S	L=150cm, 2mm safety plug (*)	
F9022D	L=60cm, safety 1,5mm DIN socket (*)	
F9052D	L=90 cm, safety 1,5 mm DIN socket (*)	
F9032D	L=120 cm, safety 1,5 mm DIN socket (*)	
F9032DN/200	L=200 cm, safety 1,5 mm DIN socket, black	
CABLES IN KIT		
ref.	description	picture
F9022D/KITX	Kit of dnap connection cables, L=60cm, safety 1,5mm DIN socket (**)	
F9052D/KITX	Kit of dnap connection cables, L=90cm, safety 1,5mm DIN socket (**)	
F9032D/KITX	Kit of dnap connection cables, L=120cm, safety 1,5mm DIN socket (***)	

(*) suffix colour key: R=red, N=black, A= Orange, V=green, G=yellow, W=white, M=brown, B=blue, I=purple

(**) F9022D/KITX and F9052D/KITX models are available with 3, 5, 7 or 8 cables, please replace "X" with the related number

(***) F9032D/KITX is available with 3, 5 or 7 cables, please replace "X" with the related number

GEL

GEL range includes models specifically intended for ECG and ultrasound: dispensers are available separately as spare parts.

ECG

ref.	description	picture
G004	spray gel, 230ml bottle	
G005	liquid gel, 260ml bottle	


ULTRASOUND

ref.	description	picture
G006	liquid gel, 260ml bottle	
G007	liquid gel, 1000ml bottle	
G0064	liquid gel, 5000ml flexible bag, 260ml refilling bottle	
G0066	liquid gel, 5000ml jerry can, 260ml refilling bottle	

ULTRASOUND AND IPL

ref.	description	picture
G002	transparent liquid gel, 20ml pouch	
G008	transparent liquid gel, 260ml bottle	
G009	transparent liquid gel, 1000ml bottle	
G010	transparent liquid gel, high density 1000ml jar	
G0084	transparent liquid gel, 5000ml flexible bag, 260ml refilling bottle	

DISPENSERS

ref.	description	picture
G0067	reusable dispenser for 5000ml jerry can	
G0068	reusable dispenser for 5000ml flexible bag	
G0069	reusable dispenser for 1000ml bottle	



FIAB SpA
Via Costoli, 4 - I-50039 Vicchio (Florence) - Italy
Tel. +39 0558497999
export@fiab.it - www.fiab.it



MD 77846